

**STAHLÉAN BRAZIL**  
SAVING LIVES



# Plastic Surgery Instruments



[www.stahleanbrazil.com](http://www.stahleanbrazil.com)



# Plastic/Oral

Scissors	3
Needle Holders	4
Forceps	5
Raspatories and Elevators	6
Raspatories	7
Spring Forceps, Hooks and Picks	9
Retractors	10
Osteotomes, Chisels and Gouges	12
Miscellaneous	13


# SCISSORS

STAHLÉAN BRAZIL  
SAVING LIVES



## Kilner Scissors

Rounded blades, 115mm long



Code	Description
Sb-PL/001a	Straight, sharp points
Sb-PL/001b	Curved, blunt points

## Aufricht Dissecting Scissors


Heavy pattern, triangular blades, 140mm long



Code	Description
Sb-PL/002a	Straight
Sb-PL/002d	Curved

## Aufricht Dissecting Scissors

Light pattern, 150mm long



Code	Description
Sb-PL/002b	Straight
SbPL/002c	Curved

## Scissors


Curved, sharp points



Code	Description
Sb-PL/004	130mm long

## Scissors


Angular, for cutting cartilage, wire cutting notch



Code	Description
Sb-PL/006	127mm long

## McIndoe Scissors


Curved on flat, round tapering ends



Code	Description
Sb-PL/007	190mm long

## Stevens Tenotomy Scissors


Supercut, 114mm long



Code	Description
Sb-OP/424	Curved/blunt
Sb-OP/425	Curved/sharp
Sb-OP/426	Straight/blunt
Sb-OP/427	Straight/sharp

## Cartilage Cutting Scissors


Fine points, serrated edges



Code	Description
SbPL/009	190mm long

## Gorney Supercut Scissors


Rhytidectomy scissor, tungsten carbide, one micro serrated and one razor sharp edge



Code	Description
Sb-PL/085a	125mm long
Sb-PL/085	200mm long, straight
Sb-PL/085b	200mm long, curved

## Fickling (Cathcart) Scissors

Angled



Code	Description
Sb-PL/049	180mm long


# NEEDLE HOLDERS

STAHLÉAN BRAZIL  
SAVING LIVES



## Foster Scissor/Needle Holder


Fenestrated jaw, 125mm long



Code	Description
Sb-PL/010a	Left

## Lawrence Needle Holder


Fine jaw, 150mm long



Code	Description
Sb-PL/012a	
Sb-PL/012b	Tungsten carbide, with two gold bows

## Derf Needle Holder

Tungsten carbide insert, with two gold bows



Code	Description
Sb-PL/013	115mm long
Sb-PL/013L	125mm long

## Reverdin Suture Needle


With automatic eye, 203mm long



Code	Description
Sb-PL/013a	Fully curved
Sb-PL/013b	Curved
Sb-PL/013c	Slightly curved
Sb-PL/013d	Straight


## Masing Needle Holder

Tungsten carbide insert, one angled bow



Code	Description
Sb-GI/009R	130mm long

## Castroviejo Needle Holder



Code	Description
Sb-OP/021	140mm long, with catch
Sb-OP/020	140mm long, without catch

### Diamond dust jaws

Code	Description
Sb-THO/060	180mm long, without catch
Sb-THO/061	200mm long, without catch
Sb-THO/062	215mm long, without catch


ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE



### Tungsten carbide jaws

Code	Description
Sb-THO/058a	180mm long
Sb-THO/058	200mm long
Sb-THO/059	215mm long

ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE




# FORCEPS

STAHLÉAN BRAZIL  
SAVING LIVES



## Mustarde Suture Holding Forceps


Convex and concave jaw



Code	Description
Sb-PL/011	133mm long

## Walsham Septum Straightening Forceps


220mm long



Code	Description
Sb-PL/030a	Straight
Sb-PL/030b	Right
Sb-PL/030c	Left

## Fickling Artery Forceps


190mm long



Code	Description
Sb-PL/044a	Serrated
Sb-PL/044b	1-2 teeth


## Bone Cutting Forceps

Straight, compound action



Code	Description
Sb-PL/033	146mm long


## Asch Septum Forceps



Code	Description
Sb-PL/031	235mm long

## Ward Bone Cutting Forceps


Slightly angled on flat



Code	Description
Sb-PL/038	178mm long

## Jansen-Middleton Septum Forceps


Compound action, 178mm long



Code	Description
Sb-PL/032	5mm bite

## Symmetry Vein Strippers

- Proven, cost e'ective removal of saphenous varicose veins
- Eliminate the need to clean and process them after use
- 30-year heritage of quality



# RASPATORIES AND ELEVATORS

## Cleft Palette Raspatory



Code	Description
Sb-PL/091	Curved downwards
Sb-PL/092	Curved upwards

## Cleft Palette Raspatory

Curved upwards



Code	Description
Sb-PL/091	Left
Sb-PL/092	Right

## McIndoe Raspatory

For cleft palette, spoon end, 170mm long



Code	Description
Sb-/PL088	Small
Sb-/PL089	Medium
Sb-/PL090	Large

## Farabeuf Elevator

Light, 160mm long



Code	Description
Sb-PL/024	12mm curved
Sb-PL/024a	8mm curved

## Howarth Elevator



Code	Description
Sb-PL/025	220mm long

## Syme Periosteal Elevator


Double ended



Code	Description
Sb-PL/046	170mm long

## Cryer Elevator

150mm long



Code	Description
Sb-PL/069	Right
Sb-PL/069a	Left

## Coupland Elevator

Curved, 150mm long



Code	Description
Sb-PL/068	3mm wide
Sb-PL/068a	3.5mm wide
Sb-PL/068b	4mm wide


# RASPATORIES

STAHLÉAN BRAZIL  
SAVING LIVES



## Aufricht Raspatory


9mm wide, 220mm long



Code	Description
Sb-PL/095	Up stroke
Sb-PL/095a	Down stroke

## Gallagher Raspatory


190mm long



Code
Sb-PL/098

## Durogrip Raspatory

7mm wide, 205mm long




12/14 teeth

Code	Description
Sb-PL/096	12/14 teeth
Sb-PL/096a	16/20 teeth
Sb-PL/096b	24/30 teeth
Sb-PL/096c	36/44 teeth
Sb-PL/096d	52/62 teeth

## Joseph Raspatory


Chequered, 165mm long



Code	Description
Sb-PL/099	Fine, 7.5mm
Sb-PL/099a	Coarse, 9mm

## Fomon Raspatory


6.5 and 7.5mm wide, 215mm long



Code	Description
Sb-PL/097	Chequered
Sb-PL/097a	Knotched

## Lewis Raspatory

8mm wide, 175mm long



Code	Description
Sb-PL/100	Fine teeth
Sb-PL/100a	Coarse teeth

# RASPATORIES


STAHLÉAN BRAZIL  
SAVING LIVES



## Maltz Raspatory

180mm long


Code	Description
Sb-PL/101	Straight, 9mm
Sb-PL/101b	Straight, 5mm
Sb-PL/101c	Curved, 5mm




## Miller Raspatory

180mm long

Code	Description
Sb-PL/102	4.3 & 6mm
Sb-PL/102a	5 & 6mm
Sb-PL/102b	5 & 8mm



## Passivation Process



### Did you know...

**Stainless steel is self passivating which provides natural corrosion resistance.**

Passivation is the natural process that occurs to a surgical instrument containing chromium. The chromium in the stainless steel naturally reacts with oxygen in the air, creating the very thin passive layer which protects the surface of the steel.

# SPRING FORCEPS, HOOKS AND PICKS

## Gillies Dissecting Forceps



Code	Description
Sb-PL/052	1 x 2 teeth, 152mm long

## McIndoe Dissecting Forceps

152mm long



Code	Description
Sb-PL/053a	1-2 teeth
Sb-PL/053b	Serrated

## Waugh Dissecting Forceps

Fine, serrated



Code	Description
Sb-GI/261	152mm long
Sb-GI/262	178mm long
Sb-GI/263	203mm long

## College Pattern Stopping and Dressing Forceps

152mm long



Code	Description
Sb-PL/048	Serrated

## Mustarde Dissecting and Suture Holding Forceps


Combined double hooked jaws



Code	Description
Sb-PL/050	156mm long

## Childe Approximating Forceps

With magazine for Michel clips



Code	Description
Sb-PL/073	180mm long

## Mitchell Trimmer

165mm long



Code
Sb-ENT/013

## Dental Pick



Code	Description
Sb-PL/067	160mm long

190mm long



Code	Description
Sb-PL/065a	Standard
Sb-PL/065b	Half Size



## Joseph Double Skin Hook

165mm long



Code
Sb-PL/056a

## Kleinert Kutz Skin Hooks



Code	Description
Sb-PL/109	16cm large
Sb-PL/108	16cm medium
Sb-PL/107	16cm small

# RETRACTORS

STAHLÉAN BRAZIL  
SAVING LIVES



## Alms Skin Retractor


Self retaining



Code	Description
Sb-PL/060	70mm long
Sb-PL/060a	100mm long

## Austin Kay Retractor

85mm long



Code	Description
Sb-PL/074	16mm wide

## Kilner Retractor


Double ended, cat's paw



Code	Description
Sb-PL/056	152mm long

## Bowdler Henry Retractor

170mm long



Code	Description
Sb-PL/076	12mm wide

## Kilner Alae Retractor


85mm long



Code	Description
Sb-PL/063a	10mm
Sb-PL/063b	13mm

## Tongue Retractor


190mm long



Code	Description
Sb-PL/071	12mm wide

## Cheek Retractor

140mm long



Code	Description
Sb-PL/078	40mm wide

## Ragnell Davis Retractor

140mm long



Code
Sb-PL/074a

symmetry surgical  
authorised distributor

**Bovie**

**Bovie ožer Electrosurgery units from 10 watts, suitable for Desiccation, all the way up to 300 units.**

We recommend the A1250S as the all round unit, for plastic surgery as it has both cut and coagulate power, up to 120 watts.




# RETRACTORS

STAHLÉAN BRAZIL  
SAVING LIVES



## Flap Retractor


180mm long



Code	Description
Sb-PL/079	32mm wide

## Kilner Cheek Retractor

Double ended, 167mm long



Code	Description
Sb-PL/058a	Small, 19mm and 25mm wide
Sb-PL/058b	Large, 29mm and 35mm wide

## Kilner Nasal Retractor


152mm long



Code	Description
Sb-PL/061a	7mm wide

## Laster Cheek Retractor


110mm long



Code	Description
Sb-PL/081	35mm wide

## Hajek Retractor


Upper lip, 165mm long



Code	Description
Sb-PL/059	20mm wide

## Freece Oral Retractor

150mm long



Code	Description
Sb-PL/080	33mm wide

## Symmetry Skin Markers

Designed to provide a permanent skin marking solution for intact skin that won't wash off while preparing a patient for a procedure

- Protects risk of cross contamination
- Sterile Skin Markers require no preoperative preparation
- Packaged with a flexible, plastic ruler that follows body contours for easy marking and reference photos
- Available in standard tip for broad clear lines as well as in micro tips for extra precision
- Packaged 36 pre-sterile Skin Markers per box

+  
Choice of  
standard tips  
or micro tips



# OSTEOTOMES, CHISELS AND GOUGES

## Ward (East Grinstead) Osteotome

Square handle, 178mm long



Code	Description
Sb-PL/040aeg	3mm wide
Sb-PL/040beg	5mm wide

140mm long



Code	Description
Sb-PL/026	3mm wide
SbPL/026a	5mm wide
Sb-PL/026b	7mm wide
Sb-PL/026c	9mm wide
Sb-PL/026g	11mm wide
Sb-PL/026d	16mm wide

## Jenkins Chisel

140mm long



Code	Description
Sb-PL/027a	2mm wide
Sb-PL/027b	4mm wide
Sb-PL/027c	6mm wide
Sb-PL/027d	8mm wide
Sb-PL/027e	10mm wide

## McIndoe Nasal Chisel

152mm long



Code	Description
Sb-PL/028a	5mm wide
Sb-PL/028b	7mm wide
Sb-PL/028c	9mm wide
Sb-PL/028d	11mm wide

## Silver Nasal Chisel

182mm long



Code	Description
Sb-PL/029a	Straight
Sb-PL/029b	Curved to right
Sb-PL/029c	Curved to left

## Pterigoideus (Maxilla) Chisel



Code	Description
Sb-PL/087	255mm long

## Tilley Septum Gouge


V shaped, straight



Code	Description
Sb-PL/043	127mm long

## Fickling Dental Mallet

191mm long




Code	Description
Sb-PL/044	25mm ø, tapering head

# MISCELLANEOUS

STAHLÉAN BRAZIL  
SAVING LIVES




## Silver Skin Graft Knife



Code	Description
Sb-PL/015a	190mm long

## Doyen Mouth Gag




Code	Description
Sb-ENT/88d	Adult
Sb-ENT/88e	Child

## Eckhož Mapping Pen



Code	Description
Sb-PL/019	145mm long
Sb-PL/019a	Nibs packs of 10

## Fergusson Ackland Mouth Gag



## Sommerland Surgical Marking Pen

Adjustable four-pointed tip




Code	Description
Sb-PL/020	152mm long

Code	Description
Sb-ENT/209	Adult, 210mm
Sb-ENT/209a	Child, 155mm


## Gillies Suction Tube

Control valve, 180mm long




Code	Description
Sb-PL/075a	2.5mm
Sb-PL/075b	3mm
Sb-PL/075c	3.5mm

## Dental Cartridge Syringe



Code	Description
Sb-PL/072	115mm long

## Watson Skin Graft Knife



Code	Description
Sb-PL/017	Right-handed
Sb-PL/017L	Left-handed

# MISCELLANEOUS

STAHLÉAN BRAZIL  
SAVING LIVES



## Freeman Areola Nipple Sizer

1.5mm

Code	Description
Sb-PL/NIP25	25mm ø
Sb-PL/NIP40	40mm ø
Sb-PL/NIP45	45mm ø
Sb-PL/NIP50	50mm ø
Sb-PL/NIP55	55mm ø



## Unna Comedo Expressor Double Ended Curette


2mm / 1mm holes

Code	Description
Sb-PL/035	135mm long



## Kelsey Fry Bone Awl

148mm long



Code	Description
Sb-PL/039a	Straight
Sb-PL/039b	Curved

## Lang Steel Scoop


Double ended

Code	Description
Sb-PL/064	7mm x 4mm, 140mm long



## Joseph Nasal Saw

Blade 35mm long



### Horizontal shaft

Code	Description
Sb-PL/034a	Cranked right
Sb-PL/034b	Cranked left

### Vertical shaft

Code	Description
Sb-PL/034c	Cranked right
Sb-PL/034d	Cranked left

## McHardy Curette

136mm long

Code	Description
Sb-PL/036	1mm, 1.5mm, and 2mm ø



## MacEwen Curette

Combined with seeker, 165mm long

Code	Description
Sb-PL/045a	Oval
Sb-PL/045b	Round



## Nero Smooth Cut Scissors

A ceramic-coated range of scissors with unbeatable durability, life cycle cost, anti-glare, cleaning efficacy, and the smoothest ever performance.





# STAHLÉAN BRAZIL

SAVING LIVES

## CONTACT US

Office in são paulo

☎ +5511972111007

☎ +5511933211007

☎ +923347711007

📷 @stahleanbrazil

✉ info@stahleanbrazil.com

🌐 www.stahleanbrazil.com