

STAHLÉAN BRAZIL
SAVING LIVES



Post-Mortem Instruments



www.stahleanbrazil.com

Post Mortem/ Autopsy

Scalpels and Knives	3
Scissors	5
Bone Cutting Forceps	7
Rib Shears	8
Dissecting Forceps	9
Forceps	10
Saws	11
Gouges	12



SCALPELS AND KNIVES

STAHLÉAN BRAZIL
SAVING LIVES



Baron Scalpel Handle


130mm long



Code	Description
Sb-PM/003	No 3
Sb-PM/004	No 4

Scalpel Handle

127mm long



Code	Description
Sb-PM/007	No 3
Sb-PM/008	No 4
Sb-PM/008a	No 4 L
Sb-PM/001a	No 7

Solid Forged Scalpel

Stainless steel
140mm long



Code	Description
Sb-PM/91	Blade length 25mm

Solid Forged Scalpel

Stainless steel
160mm long



Code	Description
Sb-PM/92	Blade length 40mm

Solid Forged Scalpel

Stainless steel
165mm long



Code	Description
Sb-PM/93	Blade length 50mm

Solid Forged Scalpel

Stainless steel
178mm long



Code	Description
Sb-PM/94	Blade length 70mm

Catlin Knife

Solid forged stainless steel



Code	Description
Sb-PM/011	Blade length 152mm

Cartilage Knife

Solid forged stainless steel



Code	Description
Sb-PM/012	Blade length 133mm
Sb-PM/013	Blade length 101mm

Resection Knife


Solid forged stainless steel



Code	Description
Sb-PM/017	Blade length 76mm

Pelvic Organ Knife

Stainless steel, curved, double edged



Code
Sb-PM/020

SCALPELS AND KNIVES

STAHLÉAN BRAZIL
SAVING LIVES



Slicing Knife

For brain or organ,
solid forged stainless steel,
22mm blade



Code	Description
Sb-PM/009NH	Blade length 254mm
Sb-PM/010NH	Blade length 152mm

Brain Knife

Solid forged stainless steel
handle, 25mm blade



Code	Description
Sb-PM/014	Blade length 254mm
Sb-PM/016	Blade length 314mm
Sb-PM/018	Blade length 254mm, 30mm blade

Organ Knife

Solid forged stainless steel,
old style, 25mm blade



Code	Description
Sb-PM/014NH	Blade length 254mm
Sb-PM/016NH	Blade length 314mm

Organ Knife and Saw

Solid forged stainless steel, 22mm blade



Code	Description
Sb-PM/017NH	Blade length 255mm

Organ Knife



Code	Description
Sb-PM/100NH	Blade length 150mm, 22mm blade


Bistoury Knife

Probe end, curved, 70mm blade



Code	Description
Sb-PM/023	176mm long

Rib Knife



Code
Sb-PM/025


SCISSORS

STAHLÉAN BRAZIL
SAVING LIVES



Scissors

Straight



Blunt/blunt

Code	Description
Sb-GI/133	130mm long
Sb-GI/136	146mm long
Sb-GI/139	165mm long


Sharp/blunt

Code	Description
Sb-GI/134	130mm long
Sb-GI/137	146mm long
Sb-GI/140	165mm long

Sharp/sharp

Code	Description
Sb-PM/034	130mm long
Sb-PM/035	146mm long
Sb-PM/036	165mm long

Dissecting Scissors



Sharp/sharp

Code	Description
Sb-PM/046	127mm long

Sharp/probe


Code	Description
Sb-PM/048	100mm long
Sb-PM/047	127mm long
Sb-PM/052	203mm long
Sb-PM/052a	180mm long

Probe/probe

Code	Description
Sb-PM/050	127mm long
Sb-PM/053	203mm long

Scissors - Mayo


Straight, round point



Code	Description
Sb-GI/409	140mm long
Sb-GI/409a	160mm long
Sb-GI/411	215mm long

Dressing Scissors


Fine end, tungsten carbide, supercut, 200mm long



Code	Description
Sb-PM/147	Blunt/sharp
Sb-PM/149	Blunt/blunt

Scissors - Mayo


Fine point



Code	Description
Sb-PM/039	Straight 139mm long
Sb-PM/043	Curved 139mm long
Sb-PM/044	Curved 190mm long
Sb-PM/045	Curved 215mm long

Artery Scissors

Small, for all arteries, tungsten carbide, supercut, 100mm long



Code	Description
Sb-PM/160	Sharp/probe


SCISSORS

STAHLÉAN BRAZIL
SAVING LIVES



Dissection Scissors


Tungsten carbide, supercut, 175mm long



Code	Description
Sb-PM/150	Blunt/blunt

Barnards Bowel Scissors


Guarded point, straight



Code	Description
Sb-PM/056	279mm long

Mayo Dissection Scissors

Tungsten carbide, supercut, for soft tissue, 175mm long



Code	Description
Sb-PM/151	Sharp/blunt

Iris scissors


102mm long



Code	Description
Sb-OP/032	Straight
Sb-OP/033	Curved

Bowel Scissors


Tungsten carbide, supercut, for dissection, straight, 180mm long



Code	Description
Sb-PM/146	Blunt/blunt

Strabismus Scissors


114mm long



Code	Description
Sb-OP/29	Straight
Sb-OP/30	Curved

Scissors


Heavy



Code	Description
Sb-PM/054	152mm long
Sb-PM/055	177mm long

Metzenbaum Scissors


185mm long



Code	Description
Sb-GI/405	Straight
SbGI/407	Curved

Bowel Scissors

Guarded point, straight



Code	Description
Sb-PM/057	196mm long


BONE CUTTING FORCEPS

STAHLÉAN BRAZIL
SAVING LIVES




Liston Bone Cutting Forceps

Straight



Code	Description
Sb-OR/032	200mm long
Sb-OR/027	220mm long
Sb-PM/066	254mm long


Ferguson "Lion" Bone Holding Forceps



Code	Description
Sb-PM/069	203mm long

Horsley Bone Cutter


Compound action
266mm long



Code	Description
Sb-PM/074	Straight
Sb-PM/075	Angled

Liston Bone Cutting Forceps

Angled



Code	Description
Sb-PM/067	254mm long
Sb-PM/068	304mm long


RIB SHEARS

STAHLÉAN BRAZIL
SAVING LIVES




Rib Shears

London Hospital pattern



Code	Description
Sb-PM/070	241mm long


Pollock Rib Shears



Code	Description
Sb-PM/071	311mm long


Rib Shears

Shark Instruments pattern




Code	Description
Sb-PM/073a	387mm long

Eslander Rib Shears



Code	Description
Sb-PM/072	215mm long

Infant Skull Shears




Code	Description
Sb-PM/073	184mm long

DISSECTING FORCEPS AND NEEDLE HOLDER

Dissecting Forceps


Serrated, turn over end



Code	Description
Sb-GI/244	178mm long

Block End Dissecting Forceps


Straight, serrated



Code	Description
Sb-GI/256	152mm long
Sb-GI/256a	203mm long
SbPM/080	254mm long
SbPM/081	310mm long

Treves Dissecting Forceps


200mm long



Code	Description
Sb-PM/084	1-2 teeth
Sb-PM/086	2-3 teeth

Jean Dissecting Forceps

177mm long




Code	Description
Sb-PM/097	Serrated
Sb-PM/098	1-2 teeth



McIndoe Dissecting Forceps

152mm long



Code	Description
Sb-GI/476	Serrated
SbGI/477	1-2 teeth

Ramsey Dissecting Forceps

177mm long

Code	Description
Sb-PM/101	Serrated
Sb-PM/102	1-2 teeth
Sb-PM/103	2-3 teeth

Mayo Hegar Needle Holder

Code	Description
Sb-GI/244	159mm long
Sb-PM/36c	178mm long
SbGI/398	203mm long


FORCEPS, CLAMPS AND RASPATORY

STAHLÉAN BRAZIL
SAVING LIVES



Spencer Wells Artery Forceps


Straight



Code	Description
Sb-GI/035	127mm long
Sb-GI/036	152mm long
Sb-GI/037	177mm long
Sb-GI/038	203mm long


Spencer Wells Artery Forceps

Curved



Code	Description
Sb-GI/039	127mm long
Sb-GI/040	152mm long
Sb-GI/041	177mm long
Sb-GI/042	203mm long


Dura Mater Stripping Forceps



Code	Description
Sb-PM/091	203mm long

Lane Bowel Clamp

Stainless steel, straight



Code	Description
Sb-PM/087	177mm long

Raspatory

For scalp



Code	Description
Sb-PM/090	175mm long

SAWS


STAHLÉAN BRAZIL
SAVING LIVES



Fixed Back Saw

14 teeth per inch

Code	Description
Sb-PM/105	Blade length 228mm
Sb-PM/107	Blade length 254mm



Barnard's Saw


14 teeth per inch

Code	Description
Sb-PM/104	Blade length 279mm

Hinged Back Saw

14 teeth per inch

Code	Description
Sb-PM/106	Blade length 228mm
Sb-PM/108	Blade length 254mm



Spinal Saw

Fixed back, double bladed, adjustable


Code	Description
Sb-PM/109	Blade length 152mm, 297mm long



Amputation Saw

Fixed back, with shaped end

Code	Description
Sb-PM/110	Blade length 254mm



GOUGES AND CHISELS

STAHLÉAN BRAZIL
SAVING LIVES




Gouge



Code	Description
Sb-PM/133	Blade 19mm

Chisel


Straight



Code	Description
Sb-PM/130	Blade 12mm
Sb-PM/131	Blade 19mm
Sb-PM/132	Blade 25mm

Spinal Chisel


With locating point



Code	Description
SbPM/135	180mm

Chisel



T shaped



Code	Description
Sb-PM/137	23mm blade, 100mm long
Sb-PM/138	25mm blade, 127mm long

Vertebrae Chisel


Angled



Code	Description
Sb-PM/136	Blade 76mm

Vertebrae/Bone Chisel

Angled, with blue acetic "Suregrip" handles



Code	Description
Sb-PM/136BH	Blade 76mm

MALLETS, PROBES, RETRACTORS, NEEDLES

Rubber Mallet


With wooden handle
305mm long



Code	Description
Sb-PM/142	300 grams

Mallet


Army pattern, lead filled



Code	Description
SbPM/141	200mm long

Hammer

Wrench end



Code	Description
Sb-PM/140	200mm long

Silver Plated Probe

With eye



Code	Description
Sb-GI/289	152cm long
Sb-GI/291	203mm long
Sb-GI/292	255mm long


Silver Probe

Fish tail



Code	Description
Sb-GI/294	152cm long
Sb-GI/295	203mm long
Sb-GI/296	255mm long


Durham Retractor



Code	Description
Sb-PM/121	215mm long

Scalpel Blade Remover

Serrated



Code	Description
Sb-PM/128	148mm long





STAHLÉAN BRAZIL

SAVING LIVES

CONTACT US

Office in são paulo

 +5511972111007

 +5511933211007

 +923347711007

 @stahleanbrazil

 info@stahleanbrazil.com

 www.stahleanbrazil.com