

STAHLÉAN BRAZIL
SAVING LIVES



Spinal Instruments



www.stahleanbrazil.com

Spinal Surgery

Forceps	3
Scissors	5
Dissectors and Probes	6
Hooks and Curettes	7
Retractors	8
Elevators	11
Suction Tubes	12
Rongeurs	13



FORCEPS



Adson Dissecting Forceps

Bayonet, 1-2 teeth, jaw 2.5mm at tip



Code	Description
SB-NEU/001c	191mm long


Jefferson Dressing Forceps



Code	Description
SB-NEU/013	Serrated jaws, 178mm long

Adson Dissecting Forceps

Bayonet, fine, serrated jaws, jaw 2.5mm at tip



Code	Description
SB-NEU/009	178mm long

Jefferson Forceps



Code	Description
SB-NEU/002	Fine, 1-2 teeth, 178mm long


Northfield Dressing Forceps

Extra light, serrated jaws



Code	Description
SB-NEU/014	1mm wide, 203mm long

Pennybacker Dissecting Forceps



Code	Description
SB-NEU/005	3-4 teeth, 165mm long

FORCEPS

STAHLÉAN BRAZIL
SAVING LIVES



Silcock Iris Forceps

102mm long



Code	Description
SB-NEU/006	1-2 teeth
SB-NEU/015	Serrated jaws

Dott Dissecting Forceps

Fine, 1-2 teeth



Code	Description
SB-NEU/004	178mm long

Cairns Dissecting Forceps


Delicate, 2-3 teeth



Code	Description
SB-NEU/003a	165mm long
SB-NEU/003b	216mm long

Dott Dressing Forceps

Angled



Code	Description
SB-NEU/018	114mm long

Cairns (Dott) Dressing Forceps


Serrated jaws



Code	Description
SB-NEU/010	Medium, 2mm wide, 216mm long
SB-NEU/012	Fine, 1mm wide, 178mm long

Tilley Dressing Forceps

Angled, serrated ends



Code	Description
SB-NEU/017	89mm long


SCISSORS

STAHLÉAN BRAZIL
SAVING LIVES



Cairns Scissors


Angled on flat



Code	Description
SB-NEU/020	178mm long

Olivecrona Scissors


Angled, with all micro serrated blade



Code	Description
SB-NEU/024	22mm blade, 150mm shaft

McKissock Scissors


Straight



Code	Description
SB-NEU/021a	146mm long

Schmieden (Taylor) Scissors

Angled, probe-pointed under blade



Code	Description
SB-NEU/023	146mm long

The Original Titanium Rhoton Dissector Set

from stahleanbrazil




+
19 dissectors
in a purpose
designed
case

- ✓ Made from lightweight titanium for balance, durability and corrosion resistance
- ✓ Design features include round, deeply knurled handles which facilitate non-slip, finger-controlled micro manipulation

Sweet Scissors

Delicate



Code	Description
SB-NEU/025a	7mm blade, 144mm shaft
SB-NEU/030	13mm blade, 117mm shaft

DISSECTORS AND PROBES

Cairns Guide and Director

Malleable



Code	Description
SB-NEU/080	127mm long

Pennybacker Probe Dissector

Double ended



Code	Description
SB-NEU/101	229mm long

Dissector

200mm long



Code	Description
SB-NEU/N30210	Soft
SB-NEU/N30220	Soft, with hole

Swedish Dissector

Double ended



Code	Description
SB-NEU/100a	Light model, 178mm long
SB-NEU/100b	Heavy model, 200mm long

Penfield Dissector



No 1

Broad curve, with 7.5mm sharp round spoon, 184mm long



No 2

7.5mm sharp round spoon with wax packer, 178mm long



No 3

Full curve and wax packer, 197mm long



No 4

Slight curve, 4mm wide, 216mm long

Code	Description
SB-NEU/069a	No 1
SB-NEU/069b	No 2
SB-NEU/069c	No 3
SB-NEU/069d	No 4

Watson Cheyne Dissector

Fine, double-ended, with probe end



Code	Description
SB-GI/287	127mm long
SB-GI/288	191mm long
SB-GI/288a	230mm long
SB-GI/288b	255mm long

HOOKS AND CURETTES

Hooks

Adson Hook


178mm long



Code	Description
SB-NEU/078a	Angled
SB-NEU/078b	Straight
SB-NEU/078b-b	Straight, blunt

Adson Oliver Neuro Hook

191mm long



Code	Description
SB-NEU/082a	Angled downward, blunt
SB-NEU/082b	Sharp
SB-NEU/082c	Angled upward, blunt

Cairns Tissue Hook

127mm long



Code	Description
SB-NEU/081b	Half round, sharp
SB-NEU/081c	Angled, sharp

Oliver Tissue Hook

Blunt



Code	Description
SB-NEU/081a	127mm long

Curettes

Falconer Scoop Curette


Double ended, angled ends



Code	Description
SB-NEU/099	191mm long

American Spinal Fusion Curettes

Blue acetal handles, 254mm long



Straight

Code	Description
SB-NEU/074a	Size 000 2mm x 4mm
SB-NEU/074b	Size 00 2.5mm x 4.5mm
SB-NEU/074c	Size 0 2.75mm x 5mm
SB-NEU/074d	Size 1 3mm x 5.5mm
SB-NEU/074e	Size 2 3.25mm x 6mm
SB-NEU/074f	Size 3 4mm x 7mm
SB-NEU/074g	Size 4 6mm x 9.5mm
SB-NEU/074h	Size 5 7mm x 14mm
SB-NEU/074i	Size 6 9mm x 17mm

Angled

Code	Description
SB-NEU/074j	Size 000 2mm x 4mm
SB-NEU/074k	Size 00 2.5mm x 4.5mm
SB-NEU/074l	Size 0 2.75mm x 5mm
SB-NEU/074m	Size 1 3mm x 5.5mm
SB-NEU/074n	Size 2 3.25mm x 6mm
SB-NEU/074o	Size 3 4mm x 7mm
SB-NEU/074p	Size 4 6mm x 9.5mm
SB-NEU/074q	Size 5 7mm x 14mm
SB-NEU/074r	Size 6 9mm x 17mm

ALSO AVAILABLE WITH ACETAL SUREGRIP HANDLE
ADD SUFFIX BH TO CATALOGUE NUMBER


RETRACTORS

STAHLÉAN BRAZIL
SAVING LIVES



Mastoid Retractor


3 x 3 blunt hooks



Code	Description
SB-NEU/N30350	114mm long

Henderson Retractor

Self retaining, straight, 4 x 4 blunt pointed blades



Code	Description
SB-NEU/030a	178mm long

Jansen (Shea) Retractor


Self retaining, blunt blades



Code	Description
SB-NEU/028	114mm long

Cairns Scalp Retractor


Curved, self retaining, swivel blades



Code	Description
SB-NEU/025b	140mm long

Cone Retractor


Self retaining, hinged arms, for use in ventriculography



Code	Description
SB-NEU/027	165mm long

Schnitker Scalp Retractor


Curved, self retaining, blunt swivel blades



Code	Description
SB-NEU/026	140mm long

Jefferson Retractor


Self retaining, curved, blunt pointed blades



Code	Description
SB-NEU/029	140mm long

Cone Laminectomy Retractor

Hinged arms, with cam rack



Code	Description
SB-NEU/051	254mm long

RETRACTORS


STAHLÉAN BRAZIL
SAVING LIVES



Cone Laminectomy Retractor

Hinged arms


Code	Description
SB-NEU/052a	Left, 165mm long
SB-NEU/052b	Right, 165mm long
SB-NEU/052c	Left, 254mm long
SB-NEU/052d	Right, 254mm long



Harvey-Jackson Laminectomy Retractor

Hinged arms, blades 38mm wide x 44mm deep


Code	Description
SB-NEU/053	305mm long



Mayo-Adson (Beckmann) Laminectomy Retractor

Hinged arms

Code	Description
SB-NEU/054	292mm long




Northfield (Cairns) Retractor

Hinged arms, one pair of blades

4 prong : 25mm wide x 16mm deep
and one pair of blades


5 prong: 32mm wide x 25mm deep

Code	Description
SB-NEU/055	292mm long



Meyerdig Hemi-Laminectomy Retractor

Code	Description
SB-NEU/216	18cm long 7" blade size 50 x 10mm
SB-NEU/219	18cm long 7" blade size 50 x 20mm
SB-NEU/217	18cm long 7" blade size 60 x 10mm
SB-NEU/220	18cm long 7" blade size 60 x 20mm
SB-NEU/218	18cm long 7" blade size 70 x 10mm
SB-NEU/221	18cm long 7" blade size 70 x 20mm
SB-NEU/224	18cm long blade left 90 x 10mm
SB-NEU/225	18cm long blades left 90 x 20mm
SB-NEU/222	18cm long blades right 90 x 10mm
SB-NEU/223	18cm long blades right 90 x 20mm




RETRACTORS

STAHLÉAN BRAZIL
SAVING LIVES



Cairns-Northfield Laminectomy Retractor


203mm arm length with three pairs of interchangeable swivelling blades, 76mm, 64mm, and 51mm



Code
SB-NEU/036

Hibbs Retractor


203mm long



Code	Description
SB-NEU/050a	Blades 10mm x 38mm
SB-NEU/050b	Blades 16mm x 54mm

Krayenbuhl Retractor


6mm wide



Code	Description
SB-NEU/5905	210mm long
SB-NEU/5905a	240mm long

Cairns 'Rake' Scalp Retractor


Sharp, 210mm long



Code	Description
SB-NEU/032a	3 prong, 16mm wide
SB-NEU/032b	4 prong, 19mm wide

Love Retractor


6mm wide



Code	Description
SB-NEU/5905b	210mm long

Cairns (Cushing) Scalp Retractor


Solid blade



Code	Description
SB-NEU/033	16mm wide, 210mm long

Love Retractor

203mm long




Code	Description
SB-NEU/077a	Straight, blade 5mm x 5mm
SB-NEU/077b	Shaft angled 90° at handle
SB-NEU/077c	Shaft angled 45° at handle

Cairns Strip Retractor

Malleable, 203mm long

8 sizes: 6, 8, 10, 11, 13, 16, 19, 25mm



Code	Description
SB-NEU/034a	20 SWG
SB-NEU/034b	18 SWG

O'Connell Retractor

Double ended



Code	Description
SB-NEU/076	216mm long


ELEVATORS

STAHLÉAN BRAZIL
SAVING LIVES



Cobb Spinal Elevator


13mm wide, knurled handle



Cobb Spinal Elevator

With blue acetal handles, 279mm long


Code	Description
SB-NEU/073a	10mm wide
SB-NEU/073b	13mm wide
SB-NEU/073c	19mm wide



ALSO AVAILABLE WITH ACETAL SUREGRIP HANDLE
ADD SUFFIX BH TO CATALOGUE NUMBER

Adson Periosteal Elevator

Serrated grip, 165mm long



straight

Code	Description
SB-OR/N028	7mm wide

curved

Code	Description
SB-OR/N029	5mm wide
SB-OR/N030	7mm wide

McKissock Spade Periosteal Elevator

Straight



Code	Description
SB-NEU/072	197mm long

Cushing Periosteal Elevator


Slightly curved end



Code	Description
SB-NEU/071	197mm long

Pennybacker Elevator

Curved, sharp, 171mm long



Code	Description
SB-OR/N025	6mm wide

SUCTION TUBES

STAHLÉAN BRAZIL
SAVING LIVES



Adson Frazier Suction Tube


Angled, fully insulated, 180mm long



Code	Description
SB-NEU/094c	8 F.G.

Retractor Suction Tube


With stilette, 220mm long



Code	Description
SB-NEU/SK001STIL	Replacement stilette

Frazier Suction Tubes


Malleable, silver, angled, 180mm long



Code	Description
SB-NEU/092a	7 F.G.
SB-NEU/092b	8 F.G.
SB-NEU/092c	9 F.G.

Dandy Rongeur


Straight



Code	Description
SB-NEU/134	3mm bite, 178mm shaft

Falconer Rongeur

Oblique jaws



Code	Description
SB-NEU/122a	Upward cutting, 3mm bite, 110mm shaft
SB-NEU/122b	Upward cutting, 5mm bite, 110mm shaft
SB-NEU/122c	Downward cutting, 3mm bite, 85mm shaft
SB-NEU/122d	Downward cutting, 5mm bite, 110mm shaft


RONGEURS

STAHLÉAN BRAZIL
SAVING LIVES



Beck Rongeur


Straight, for intervertebral disc



Code	Description
SB-NEU/125	2mm bite, 130mm shaft

Love-Grunwald Rongeur


Cup shaped, jaws 3mm x 10mm, 178mm shaft



Code	Description
SB-NEU/131a	Straight
SB-NEU/131b	Curved upward
SB-NEU/131c	Curved downward

Cushing Rongeur


130mm shaft



Code	Description
SB-NEU/127-2	2mm bite
SB-NEU/127-3	3mm bite
SB-NEU/127-4	4mm bite
SB-NEU/127-5	5mm bite
SB-NEU/127-6	6mm bite

Mackay Rongeur


140°, straight, 114mm shaft



Code	Description
SB-SH/4249	5.5mm

Cushing Rongeur

130mm shaft



150° upward cutting


Code	Description
SB-NEU/128c	2mm bite
SB-NEU/128	3mm bite
SB-NEU/128a	4mm bite
SB-NEU/128b	5mm bite

150° downward cutting

Code	Description
SB-NEU/130c	2mm bite
SB-NEU/130	3mm bite
SB-NEU/130a	4mm bite
SB-NEU/130b	5mm bite

Northfield Rongeur

Curved, upward cutting



Code	Description
SB-NEU/132	4mm bite, 133mm shaft


RONGEURS

STAHLÉAN BRAZIL
SAVING LIVES



Pennybacker Rongeur


120mm shaft



Code	Description
SB-NEU/133	3mm bite
SB-NEU/133-4	4mm bite

Daniel Rongeur


Miniature, for neuropathology



Code	Description
SB-NEU/104	1.5mm bite, 127mm long

St Bart's Rongeur


114mm shaft



Code	Description
SB-ENT/140a	3mm bite

Cairns Rongeur


Fine jaws, angled on flat, serrated jaws, curved handles



Code	Description
SB-NEU/103	4mm bite 152mm long

St Bart's Rongeur


140°, downward cutting, 114mm shaft



Code	Description
SB-ENT/140ac	3mm

Danish Pattern Rongeur

Angled on flat, fine jaws, curved handles, screw joint



Code	Description
SB-NEU/109	4mm bite, 203mm long

St Bart's Rongeur


140°, upward cutting, 114mm shaft



Code	Description
SB-ENT/140b	3mm

Falconer Rongeur

Angled to side, light, screw joint



Code	Description
SB-NEU/110	6mm bite, 203mm long


RONGEURS

STAHLÉAN BRAZIL
SAVING LIVES



Wilms Rongeur


Angled on flat, narrow jaws



Code	Description
SB-NEU/114	4mm bite, 203mm long

Pennybacker Rongeur


Straight, fine jaws, 203mm long



Code	Description
SB-NEU/107	4mm bite

Glasgow Pattern Rongeur


Slightly curved on flat, screw joint



Code	Description
SB-NEU/111	4mm bite, 146mm long

Northfield Rongeur


Straight, 10mm bite



Code	Description
SB-NEU/106	One shallow blade, 200mm long

Heanley Rongeur


Long fine straight jaws, screw joint



Code	Description
SB-NEU/112	4mm bite, 178mm long

Stookey Rongeur


Angled to side, fine jaws



Code	Description
SB-NEU/108	3mm bite, 178mm long

Jansen Rongeur


Short curve, screw joint



Code	Description
SB-NEU/105	8mm bite, 152mm long

Trotter Rongeur

Angled to side, strong, screw joint



Code	Description
SB-NEU/113	10mm bite, 203mm long



RONGEURS

STAHLÉAN BRAZIL
SAVING LIVES



Luer-Jansen Rongeur


Light, curved on flat, compound action



Code	Description
SB-NEU/116	4mm bite, 197mm long

Olivecrona Rongeur


Angled to side, compound action, 229mm long



Code	Description
SB-NEU/119a	5mm bite
SB-NEU/119b	6mm bite

Ruskin Rongeur

Curved on flat, compound action



Code	Description
SB-NEU/118	6mm bite, 184mm long



RONGEURS

STAHLÉAN BRAZIL
SAVING LIVES



Leksell Rongeur



Curved to side, compound action



Code	Description
SB-NEU/115a	Small, 6mm bite, 229mm long
SB-NEU/115b	Large, 8mm bite, 254mm long

Luer Stille Bone Rongeur


Curved, 250mm long



Code	Description
SB-OR/071b	11mm bite
SB-NEU/117	Heavy, 229mm long

Swedish Pattern Laminectomy Shears

Double crank, compound action



Code	Description
SB-NEU/121	267mm long

RONGEURS

STAHLÉAN BRAZIL
SAVING LIVES

Kerrison Rongeurs

Cloward disc punch (not illustrated)
130° upward cutting, 190mm long


BLACK CERAMIC COATED KERRISONS AVAILABLE
PLEASE ADD SUFFIX -BK WHEN ORDERING

Code	Description
SB-NEU/CL002	2mm wide
SB-NEU/CL003	3mm wide
SB-NEU/CL004	4mm wide
SB-NEU/CL005	5mm wide
SB-NEU/CL006	6mm wide



Kerrison punch

1mm 2mm 3mm 4mm 5mm



Non-Detachable

90° upward cutting

Code	Description
SB-NEU/136	180mm(L) x 1mm(W)
SB-NEU/137	180mm(L) x 2mm(W)
SB-NEU/138	180mm(L) x 3mm(W)
SB-NEU/139	180mm(L) x 4mm(W)
SB-NEU/140	180mm(L) x 5mm(W)
SB-NEU/141	230mm(L) x 2mm(W)
SB-NEU/142	230mm(L) x 3mm(W)
SB-NEU/143	230mm(L) x 4mm(W)
SB-NEU/144	230mm(L) x 5mm(W)

90° downward cutting

Code	Description
SB-NEU/154	180mm(L) x 1mm(W)
SB-NEU/155	180mm(L) x 2mm(W)
SB-NEU/156	180mm(L) x 3mm(W)
SB-NEU/157	180mm(L) x 4mm(W)
SB-NEU/158	180mm(L) x 5mm(W)
SB-NEU/159	230mm(L) x 2mm(W)
SB-NEU/160	230mm(L) x 3mm(W)
SB-NEU/161	230mm(L) x 4mm(W)
SB-NEU/162	230mm(L) x 5mm(W)

130° upward cutting

Code	Description
SB-NEU/172	180mm(L) x 1mm(W)
SB-NEU/173	180mm(L) x 2mm(W)
SB-NEU/174	180mm(L) x 3mm(W)
SB-NEU/175	180mm(L) x 4mm(W)
SB-NEU/176	180mm(L) x 5mm(W)
SB-NEU/177	200mm(L) x 3mm(W)
SB-NEU/178	200mm(L) x 4mm(W)
SB-NEU/179	200mm(L) x 5mm(W)
SB-NEU/180	200mm(L) x 6mm(W)
SB-NEU/181	230mm(L) x 2mm(W)
SB-NEU/182	230mm(L) x 3mm(W)
SB-NEU/183	230mm(L) x 4mm(W)
SB-NEU/184	230mm(L) x 5mm(W)

130° downward cutting

Code	Description
SB-NEU/198	180mm(L) x 1mm(W)
SB-NEU/199	180mm(L) x 2mm(W)
SB-NEU/200	180mm(L) x 3mm(W)
SB-NEU/201	180mm(L) x 4mm(W)
SB-NEU/202	180mm(L) x 5mm(W)
SB-NEU/203	230mm(L) x 2mm(W)
SB-NEU/204	230mm(L) x 3mm(W)
SB-NEU/205	230mm(L) x 4mm(W)
SB-NEU/206	230mm(L) x 5mm(W)

Detachable

90° upward cutting

Code	Description
SB-NEU/145	180(L) x 1mm(W)
SB-NEU/146	180(L) x 2mm(W)
SB-NEU/147	180(L) x 3mm(W)
SB-NEU/148	180(L) x 4mm(W)
SB-NEU/149	180(L) x 5mm(W)
SB-NEU/150	230(L) x 2mm(W)
SB-NEU/151	230(L) x 3mm(W)
SB-NEU/152	230(L) x 4mm(W)
SB-NEU/153	230(L) x 5mm(W)

90° downward cutting

Code	Description
SB-NEU/163	180(L) x 1mm(W)
SB-NEU/164	180(L) x 2mm(W)
SB-NEU/165	180(L) x 3mm(W)
SB-NEU/166	180(L) x 4mm(W)
SB-NEU/167	180(L) x 5mm(W)
SB-NEU/168	230(L) x 2mm(W)
SB-NEU/169	230(L) x 3mm(W)
SB-NEU/170	230(L) x 4mm(W)
SB-NEU/171	230(L) x 5mm(W)

130° upward cutting

Code	Description
SB-NEU/185	180(L) x 1mm(W)
SB-NEU/186	180(L) x 2mm(W)
SB-NEU/187	180(L) x 3mm(W)
SB-NEU/188	180(L) x 4mm(W)
SB-NEU/189	180(L) x 5mm(W)
SB-NEU/190	200(L) x 3mm(W)
SB-NEU/191	200(L) x 4mm(W)
SB-NEU/192	200(L) x 5mm(W)
SB-NEU/193	200(L) x 6mm(W)
SB-NEU/194	230(L) x 2mm(W)
SB-NEU/195	230(L) x 3mm(W)
SB-NEU/196	230(L) x 4mm(W)
SB-NEU/197	230(L) x 5mm(W)

130° downward cutting

Code	Description
SB-NEU/207	180(L) x 1mm(W)
SB-NEU/208	180(L) x 2mm(W)
SB-NEU/209	180(L) x 3mm(W)
SB-NEU/210	180(L) x 4mm(W)
SB-NEU/211	180(L) x 5mm(W)
SB-NEU/212	230(L) x 2mm(W)
SB-NEU/213	230(L) x 3mm(W)
SB-NEU/214	230(L) x 4mm(W)
SB-NEU/215	230(L) x 5mm(W)

Stahleanbrazil

Your biggest challenges
solved with one product:

- Productive bite every time
- Never have to use a dull Kerrison again
- Minimise the risk of dural tears
- Reduce hand and arm fatigue
- Easily opened for cleaning
- No need to open multiple trays





STAHLÉAN BRAZIL


SAVING LIVES

CONTACT US

Office in são paulo

 +5511972111007

 +5511933211007

 +923347711007

 @stahleanbrazil

 info@stahleanbrazil.com

 www.stahleanbrazil.com