

STAHLÉAN BRAZIL
SAVING LIVES



**Cardiovascular Instruments
&
Thoracic Instruments**



www.stahleanbrazil.com

Cardiovascular Instruments & Thoracic Instruments

Artery Clamps	3
Clamps	5
Bulldog Clamps	10
Scissors	11
Needle Holders	12
Retractors	13
Rib Spreaders	14
Rib Shears	17
Elevators	18



ARTERY CLAMPS



An Overview On

Atraumatic Jaws

Also known as non-traumatic, atraumatic jaws will be mainly found in delicate surgery areas such as vascular and cardiothoracic. The design is such, that the surface area gripped is reduced, ensuring less trauma without the loss of grip.

Atraumatic jaws are available as straight or curved.

Atraumatic Jaw, Straight




Atraumatic Jaw, Curved



Cooley Artery Clamps

Box joint, atraumatic, curved



Code	Description
Sb-GI/071	160mm long
Sb-GI/072	200mm long

Large, curved 60°

Code	Description
Sb-GI/073	150mm long

Large, strongly curved 90°


Code	Description
Sb-GI/074	140mm long

Curved 30°

Code	Description
Sb-GI/075	160mm long

DeBakey Bainbridge Clamps

Box joint, atraumatic, haemostatic




180mm long

Code	Description
Sb-GI/078	Straight
Sb-GI/077	Curved

DeBakey Bainbridge Clamps

Box joint, atraumatic




Curved

Code	Description
Sb-GI/079	165mm long

DeBakey Bainbridge Clamps

Box joint, atraumatic, haemostatic




155mm long

Code	Description
Sb-GI/105	Curved
Sb-GI/106	Straight

DeBakey Clamps

Box joint, atraumatic



Curved


Code	Description
Sb-GI/086	230mm long

ARTERY CLAMPS



DeBakey Clamps

Box joint, atraumatic




Curved

Code	Description
Sb-GI/088	230mm long
Sb-GI/089	240mm long
Sb-GI/090	245mm long

DeBakey Clamps

Box joint, atraumatic




Curved

Code	Description
Sb-GI/094	150mm long

DeBakey Clamps

Box joint, atraumatic




Curved

Code	Description
Sb-GI/091	250mm long
Sb-GI/092	270mm long
Sb-GI/093	280mm long

DeBakey Clamps

Box joint, atraumatic




Curved

Code	Description
Sb-GI/096	180mm long
Sb-GI/097	203mm long

DeBakey Clamps


Box joint, atraumatic, straight



Curved

Code	Description
Sb-GI/100	160mm long
Sb-GI/104	165mm long
Sb-GI/101	203mm long
Sb-GI/098	215mm long
Sb-GI/102	220mm long
Sb-GI/103	230mm long
Sb-GI/099	245mm long

Box joint, atraumatic




Curved

Code	Description
Sb-GI/080	165mm long
Sb-GI/081	190mm long
Sb-GI/082	203mm long
Sb-GI/083	215mm long
Sb-GI/084	220mm long

CLAMPS




Duval Lung Grasping Clamps



Code	Description
Sb-THO/033	185mm long


Tudor Edwards Artery Clamps

Curved, serrated jaws, screw joint, satin finish



Code	Description
Sb-THO/015	230mm long


Babcock Lung Grasping Clamps



Code	Description
Sb-THO/034	245mm long

O'Shaughnessy Artery Clamps

Right angled serrated jaw



Code	Description
SbTHO/036a	150mm long
Sb-THO/036b	180mm long
Sb-THO/036c	200mm long
Sb-THO/036d	230mm long

**Stahleanbrazil now
offer a VATS range
of instruments
through Reuchlen**




CLAMPS



DeBakey Atraumatic Dissecting Clamps

Straight, atraumatic jaws



1.5mm

Code	Description
Sb-GI/238a	150mm long
Sb-GI/240a	240mm long

2mm

Code	Description
SI-GI/238	150mm long
SI-GI/241	180mm long

2.5mm


Code	Description
Sb-GI/238b	150mm long
Sb-GI/240b	240mm long

2.8mm

Code	Description
SI-GI/238c	150mm long
SI-GI/240c	240mm long
SI-GI/241c	180mm long

Magill Intratracheal Catheter Clamps



Screw joint




Code	Description
Sb-THO/032	Infant, 160mm long
Sb-THO/032a	Child, 190mm long
Sb-THO/032b	Adult, 230mm long

SEMB Ligature Passing Clamps

240mm long




Code	Description
Sb-THO/044	Slightly curved
Sb-THO/045	Fully curved



Lambert-Kay Clamp


Atraumatic



Code	Description
Sb-THO/029a	Curved, 220mm long
Sb-THO/029b	Straight, 215mm long
Sb-THO/029c	Straight, 200mm long

Leland-Jones Clamp

Atraumatic, hemostatic




Code	Description
Sb-THO/023a	Curved, 190mm long
Sb-THO/023b	Straight, 195mm long

CLAMPS

STAHLÉAN BRAZIL
SAVING LIVES

DeBakey Clamp


Straight jaws, ring handle



Code	Description
Sb-THO/048	137mm long

DeBakey Clamp


'S' shaped curved jaws, curved shaft, ring handle



Code	Description
Sb-THO/049	127mm long


DeBakey "S" Shaped Clamp

300mm long



Code	Description
Sb-THO/096	115mm jaws

DeBakey Atraumatic Clamp



Code	Description
Sb-THO/090	110mm jaw, 300mm long
Sb-THO/092	95mm jaw, 250mm long


CLAMPS

STAHLÉAN BRAZIL
SAVING LIVES



DeBakey Atraumatic Clamp


Right-angled, 40mm jaw



Code	Description
Sb-THO/044a	220mm long

DeBakey Atraumatic Clamp


Right-angled, 50mm jaw



Code	Description
Sb-THO/045	200mm long

DeBakey Atraumatic Clamp


Right-angled, 65mm jaw



Code	Description
Sb-THO/046	215mm long

DeBakey Atraumatic Clamp

Right-angled, 100mm jaw



Code	Description
Sb-THO/047	280mm long

CLAMPS



STAHLÉAN BRAZIL
SAVING LIVES

DeBakey Coarctation Clamp

Angled, 225mm long

Code

Sb-THO/052



DeBakey Multipurpose Clamp



Code


Sb-THO/095

Description

85mm jaw, 230mm long

Sb-THO/064

100mm jaw, 280mm long





DeBakey Multipurpose Clamp


Obtuse angled 60°, 280mm long

Code

Sb-THO/053



DeBakey 45° Paediatric Clamp




Code

Sb-THO/050

Description

130mm long

DeBakey 90° Paediatric Clamp



Code

Sb-THO/051

Description

105mm long

BULLDOG CLAMPS



Dieffenbach Bulldog Clamp



Straight, serrated

Code	Description
Sb-GI/112	1.2mm fine jaw, 40mm long
Sb-GI/113	1.2mm fine jaw, 45mm long
Sb-GI/108	1.8mm jaw, 30mm long
Sb-GI/109	1.8mm jaw, 35mm long
Sb-GI/110	2.5mm jaw, 40mm long
Sb-GI/111	2.5mm jaw, 45mm long
Sb-GI/107	6mm heavy jaw, 60mm long
Sb-GI/455	20mm jaw, 75mm long
Sb-GI/455a	30mm jaw, 85mm long
Sb-GI/455b	45mm jaw, 105mm long
Sb-GI/455c	65mm jaw, 120mm long

Straight, plain

Code	Description
Sb-GI/129	1.2mm fine jaw, 40mm long
Sb-GI/131	1.2mm fine jaw, 45mm long
Sb-GI/121	1.8mm jaw, 30mm long
Sb-GI/123	1.8mm jaw, 35mm long
Sb-GI/125	2.5mm jaw, 40mm long
Sb-GI/127	2.5mm jaw, 45mm long


Curved, serrated

Code	Description
Sb-GI/119	1.2mm fine jaw, 40mm long
Sb-GI/120	1.2mm fine jaw, 45mm long
Sb-GI/115	1.8mm jaw, 30mm long
Sb-GI/116	1.8mm jaw, 35mm long
Sb-GI/117	2.5mm jaw, 40mm long
Sb-GI/118	2.5mm jaw, 45mm long
Sb-GI/114	6mm heavy jaw, 60mm long
Sb-GI/456	20mm jaw, 75mm long
Sb-GI/456a	30mm jaw, 85mm long
Sb-GI/456b	45mm jaw, 105mm long
Sb-GI/456c	65mm jaw, 120mm long

Curved, plain

Code	Description
Sb-GI/130	1.2mm fine jaw, 40mm long
Sb-GI/132	1.2mm fine jaw, 45mm long
Sb-GI/122	1.8mm jaw, 30mm long
Sb-GI/124	1.8mm jaw, 35mm long
Sb-GI/126	2.5mm jaw, 40mm long
Sb-GI/128	2.5mm jaw, 45mm long

Blalock Bulldog Clamp



Serrated, with hole in handle

Straight

Code	Description
Sb-GI/480	11mm jaw, 38mm long
Sb-GI/482	14mm jaw, 48mm long
Sb-GI/493	17mm jaw, 60mm long
Sb-GI/481	27mm jaw, 70mm long

Curved

Code	Description
Sb-GI/478	11mm jaw, 38mm long
Sb-GI/479	14mm jaw, 48mm long
Sb-GI/491	17mm jaw, 60mm long
Sb-GI/492	27mm jaw, 70mm long

Plain, with hole in handle

Straight


Code	Description
Sb-GI/486	11mm jaw, 38mm long
Sb-GI/487	14mm jaw, 48mm long
Sb-GI/488	17mm jaw, 60mm long
Sb-GI/489	27mm jaw, 70mm long

Curved

Code	Description
Sb-GI/483	11mm jaw, 38mm long
Sb-GI/484	14mm jaw, 48mm long
Sb-GI/485	17mm jaw, 60mm long
Sb-GI/490	27mm jaw, 70mm long

Diethrich Bulldog Clamp

Angled 12mm jaw




Code	Description
Sb-JD-GI/8097	44mm long

SCISSORS

STAHLÉAN BRAZIL
SAVING LIVES

Potts Scissors


Satin finish, 190mm long



Code	Description
Sb-THO/019	Angled on flat 45°
Sb-THO/020a	Angled to side 45°
Sb-THO/020b	Acute angle to side

Potts Diethrich Scissors

Supercut, 190mm long



Code	Description
Sb-THO/078	25° angled on side
Sb-THO/079	45° angled on side
Sb-THO/080	90° angled on side

Diethrich Circumflex Artery Scissors


Supercut, 175mm long



Code	Description
Sb-THO/077	125° angled on side

Boyd Scissor


Supercut, slight curve



Code	Description
Sb-THO/047	180mm long

Lobectomy Scissors


Curved on flat



Code	Description
Sb-THO/001b	305mm long

Nelson Lobectomy Scissor

Tungsten carbide, curved on flat, 230mm long




Code
Sb-THO/057

NEEDLE HOLDERS



Needle Holder and Wire Twister

For sternum wire, left-handed, 300mm long




Code

Sb-THO/083

Julian Needle Holder

Light, delicate model, with tungsten carbide jaws
220mm long



Code

Sb-THO/081


Sb-THO/081a

Description

Left-handed

Rochester Needle Holder

Light, delicate model, with tungsten carbide jaws
190mm long



Code


Sb-THO/082

Sb-THO/082a

Description


Left-handed

DeBakey Needle Holder



Stratte Needle Holder

210mm long



Code

THO/SI-5642b

Code

THO/SI-5642a

Sb-GI/502A

Sb-GI/502

SI003944

Sb-5642A

Sb-GI/036BM

Sb-GI/037BM

Sb-GI/038BM

Description

250mm long

165mm very delicate

180mm tungsten carbide

260mm long straight with tungsten carbide inserts

240mm long


Tungsten carbide delicate 230mm long

Tungsten carbide delicate 250mm long

Tungsten carbide delicate 305mm long

Castroviejo Needle Holder

Code	Description
Sb-OP/021	140mm long, with catch
Sb-OP/020	140mm long, without catch



Diamond dust jaws


Code	Description
Sb-THO/060	180mm long, without catch
Sb-THO/061	200mm long, without catch
Sb-THO/062	215mm long, without catch



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Tungsten carbide jaws

Code	Description
Sb-THO/058a	180mm long
Sb-THO/058	200mm long
Sb-THO/059	215mm long



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

RETRACTORS



Blanco Valve Leaflet Retractor

280mm long

Code

Sb-THO/041



Harrington Splanchnic Retractor


Small size, end of blade
38mm wide x 127mm deep
satin finish handle

Code

Sb-THO/011


Description

320mm long



Brocks Erecton Spinal Back End Retractor

215mm long



Code

Sb-THO/054

Description

13 x 28mm

Cleland Valve Leaflet Retractor

250mm long



Code

Sb-THO/097

Ross Retractor

250mm long



Code

Sb-THO/089

Description

16mm x 15mm blade

Sb-THO/089a

23mm x 13mm blade

Sb-THO/089b

24mm x 16mm blade

Code

Sb-THO/089c

Description

27mm x 13mm blade

Sb-THO/089d


27mm x 27mm blade

Sb-THO/089e

25mm x 38mm blade

Tudor Edwards Retractor

25mm Wide



Code


Sb-THO/013

Description

205mm long

Allison Lung Retractor

Satin finish handle



Code

Sb-THO/012a

Description

Infant

Sb-THO/012b


Child, 280mm blade

Sb-THO/012c

Adult, 300mm blade

Cooley Atrial Retractor

260mm overall length



Code

Sb-THO/046

Description

Right

Sb-THO/046a

Left


RIB SPREADERS

STAHLÉAN BRAZIL
SAVING LIVES



Finochietto Rib Spreader


Straight arms, 20mm square blades



Code	Description
Sb-THO/006	Infant, maximum opening 80mm x 60mm long
Sb-THO/037	Adult, maximum opening 200mm, blades 44 x 65mm

Morse-Favaloro Rib Spreader

Curved arms, 275mm long



Code	Description
Sb-THO/103	Swivel blades 215mm long
Sb-THO/103a	Swivel blades 185mm long

IMA Rib Spreader

Converts left or right


Code	Description
Sb-THO/102	Standard frame only
Sb-THO/102a	Ministernotomy frame only
Sb-THO/102b	Complete, standard, 1 x frame, 1 x 100mm long x 20mm wide U shaped blade, 1 x claw blades 28mm and 32mm deep
Sb-THO/102c	Complete, ministernotomy, 1 x frame, 1 x 30mm long x 20mm wide U shaped blade, 1 x claw blades 43mm and 50mm deep

Claw blade only

Code	Description
Sb-THO/102d	28mm deep
SbTHO/102e	32mm deep
SbTHO/102f	38mm deep
Sb-THO/102g	43mm deep
SbTHO/102h	50mm deep


U shaped blade only

Code	Description
SbTHO/102i	100mm long x 20mm wide
SbTHO/102j	100mm long x 25mm wide
SbTHO/102k	30mm long x 20mm wide



Price Thomas Rib Spreader


Code	Description
Sb-THO/99	Complete with 3 pairs of blades, key, 250, 215, 68mm
Sb-THO/99a	Blade only, 50mm x 38mm, pair
Sb-THO/99b	Blade only, 50mm x 45mm, pair
Sb-THO/99c	Blade only, 50mm x 65mm, pair
Sb-THO/99d	Nut and bolt for blade only
SbTHO/99e	Frame only
Sb-THO/99f	Key only




RIB SPREADERS



Sellors Rib Spreader



Code	Description
Sb-THO/100	Complete with 3 pairs of blades, key,
Sb-THO/100a	Blade only, 50mm x 38mm, pair
Sb-THO/100b	Blade only, 50mm x 45mm, pair
Sb-THO/100c	Blade only, 50mm x 65mm, pair
Sb-THO/100d	Frame only
Sb-THO/100e	Nut and bolt for blade only
Sb-THO/100f	Key only



Harefield sternal blade only, pair


Code	Description
Sb-THO/100g	70 x 30mm
Sb-THO/100h	120 x 30mm

Harefield scapula blade only, pair

Code	Description
Sb-THO/100i	75 x 70mm

Ronald Edwards Rib Spreader


Curved arms, fenestrated



Code	Description
Sb-THO/104	Infant, 100mm, 25mm x 30mm blades
SbSH-5158	Child, 150mm, 30mm x 50mm blades
Sb-SH-5720	Adult, 245mm, 65mm x 75mm blades

Bailey Rib Approximator and Spreader


Combined, satin finish



Code	Description
Sb-THO/003	160mm long

Sellors Rib Approximator and Spreader

Combined action, satin finish




Code	Description
Sb-THO/002	195mm long

RIB SPREADERS




Cooley Rib Spreader




Code	Description
Sb-43950-50	Adult, 40mm x 110mm blades, complete with cosgrove mounting bars
Sb-43950-90	Adult, 40mm x 110mm blades, complete with cosgrove mitral valve retractor attachment, 5 x blades, 3 x clamp, 1 x side arm

Wells Rib Spreader




Code	Description
Sb-44080-65	Small rake retractor on shaft
Sb-44080-80	Frame only, 90mm wide x 25mm deep blades
Sb-44080-90	Complete, frame, pillar clamp, 4 x ball clamps, 2 x mounting bars, hex keys, 1 x tricuspid rake ret, 2 x small, medium and large rake ret

Tudor Edwards Rib Spreader




Code	Description
Sb-THO/101	Complete, 2 of each blade, plus key
Sb-THO/101a	25mm blade only, pair
Sb-THO/101b	50mm blade only, pair
Sb-THO/101c	75mm blade only, pair

Tußer Rib Spreader



Code	Description
Sb-THO/005	Blades 25mm deep x 19mm wide


Satin finish handles



Code	Description
Sb-THO/010	240mm long

RIB SHEARS AND LUNG FORCEPS


STAHLÉAN BRAZIL
SAVING LIVES



Code	Description
Sb-THO/041	270mm long

Tudor Edwards Rib Shears


Compound action, anterior



Code	Description
Sb-THO/042	235mm long

Mayo Lung Forceps

Fenestrated, semi-circular jaw, satin finish



Code	Description
Sb-THO/014	230mm long

Special Care Needed For Your Thoracic Instruments

Staining



Beware of the use of Saline during surgery. This can be very damaging to any instrument if not immediately rinsed after surgery. Saline is highly corrosive. The image shows damage to the passive surface of an instrument.

Rusting




Due to the nature of their use, Thoracic Retractors can tarnish or present rust easily. These retractors can be fully refurbished and passivated by stahleanbrazil

ELEVATORS AND MISCELLANEOUS



Tudor Edwards Periosteal Elevator


First rib, 175mm long



Code	Description
Sb-THO/009b	18mm wide

Doyen Periosteal Elevator

Satin finish handle



165mm long

Code	Description
Sb-OR/N031	Right
Sb-OR/N032	Left

180mm long

Code	Description
SbOR/256	Right
SbOR/255	Left

Blalock Hook

249mm long



Code	Description
Sb-THO/042	45° angle, 5mm tip
SbTHO/043b	90° angle, 5mm tip

Ross (Papworth) Ventricular Sucker


33mm long



Code	Description
SbTHO/PAP001	Bullet

Ross (Papworth) Sump Sucker


42mm long



Code	Description
Sb-THO/PAP002	Round end

Ross (Papworth) Flat Sucker


42mm long



Code
SI-THO/PAP003


Ross Probe

200mm long



Code	Description
Sb-THO/055	1mm
Sb-THO/055a	1.5mm

Tibbs Arterial Cannula



GU pattern

Code	Description
Sb-TIB/001	Set of four
Sb-TIB/002	Large cone
Sb-TIB/003	Large olive
Sb-TIB/004	Small cone
Sb-TIB/005	Small olive

Code	Description
Sb-TIB/006	Large cone
Sb-TIB/007	Large olive
Sb-TIB/008	Small cone
Sb-TIB/009	Small olive
Sb-TIB/010	Set of four



STAHLÉAN BRAZIL

SAVING LIVES

CONTACT US

Office in são paulo

☎ +5511 972111007

📞 +92-3347711007

📷 @stahleanbrazil

✉ info@stahleanbrazil.com

🌐 www.stahleanbrazil.com