

STAHLÉAN BRAZIL


SAVING LIVES



ENT Instruments



www.stahleanbrazil.com



ENT

Middle Ear Instruments	3
Picks and Scoops	5
Hooks and Probes	6
Retractors	7
Suction Tubes	8
Speculum	9
Aural Forceps	11
Artery Forceps	13
Tonsil, Turbinate and Septum Forceps	14
Dressing Forceps	15
Other Forceps	16
Mouth Gags	17
Tongue Plates and Depressors	18
Punches	19
Rongeurs	23
Scissors	25
Gouges, Knives and Elevators	27
Curettes and Myringotomes	28
Elevators and Dissectors	29
Miscellaneous	30




MIDDLE EAR INSTRUMENTS



Ironside Instruments


160mm long, black finish



Code	Description
Sb-ENT/200a	Knife 1, vertical incision
Sb-ENT/200b	Knife 2, horizontal incision
Sb-ENT/200c	Elevator
Sb-ENT/200d	Needle point
Sb-ENT/200e	Hook, superior elevation, angle 45°
Sb-ENT/200f	Hook, inferior elevation, angle 45°
Sb-ENT/200g	Hook, anterior elevation, angle 90° right
Sb-ENT/200h	Hook, anterior elevation, angle 90° left
Sb-ENT/200i	Measuring device

Rosen Instruments


160mm long, black finish



Code	Description
Sb-ENT/202a	Curette, tympanic ring, small
Sb-ENT/202b	Curette, tympanic ring, large
Sb-ENT/202c	Probe, footplate movement test, fine
Sb-ENT/202d	Mobiliser, slight curve
Sb-ENT/202e	Mobiliser, full curve
Sb-ENT/202f	Scoop, size 000
Sb-ENT/202g	Scoop, size 00
Sb-ENT/202h	Scoop, size 0
Sb-ENT/202i	Scoop, size 1
Sb-ENT/202j	Fenestrator, oval window

Schuknecht Instruments


160mm long, black finish



Code	Description
Sb-ENT/201a	Hook upward, 0.4mm
Sb-ENT/201b	Hook downward, 0.4mm
Sb-ENT/201c	Hook upward, 0.7mm
Sb-ENT/201d	Hook downward, 0.7mm
Sb-ENT/201e	Elevator, oval
Sb-ENT/201f	Hook, right-angled upward
Sb-ENT/201g	Pick, straight
Sb-ENT/201h	Spatula

Beale Instruments

165mm long, satin finish




Code	Description
Sb-ENT/204a	Tympanoplasty probe, fine blunt end
Sb-ENT/204b	Raspatory, slight curve
Sb-ENT/204c	Raspatory, slight curve, small
Sb-ENT/204d	Rugine, double ended, sharp
Sb-ENT/204e	Rugine, double ended
Sb-ENT/204f	Spatula, graft smoothing

MIDDLE EAR INSTRUMENTS



Beale Instruments


170mm long, satin finish



Code	Description
Sb-ENT/203a	Knife, Stapedius tendon cutting
Sb-ENT/203b	Raspatory, mucous removal
Sb-ENT/203c	Osteotome, footplate
Sb-ENT/203d	Perforator, elevation fixed footplate
Sb-ENT/203e	Hook, removal of footplate, angled downward
Sb-ENT/203f	Hook, removal of footplate, angled upward
Sb-ENT/203g	Hook, bone removal from oval window, angled left
Sb-ENT/203h	Applicator, vein graft for right ear
Sb-ENT/203i	Applicator, vein graft for left ear
Sb-ENT/203j	Curette, removal bony rim
Sb-ENT/203k	Curette

Hough Instruments


165mm long, stainless steel



Code	Description
Sb-ENT/205a	Excavator, "whirlybird", right, large
Sb-ENT/205b	Excavator, "whirlybird", left, large
Sb-ENT/205c	Excavator, "whirlybird", right, small
Sb-ENT/205d	Excavator, "whirlybird", left, small
Sb-ENT/205e	Pick, stapedectomy footplate, angled 90°, 0.33mm
Sb-ENT/205f	Pick, stapedectomy footplate, angled 45°, 0.5mm
Sb-ENT/205g	Pick, stapedectomy footplate, angled 15°, 0.5mm
Sb-ENT/205h	Hook, "hoe", excavator, angled 90°
Sb-ENT/205i	Hook, "hoe", excavator, angled 45°
Sb-ENT/205j	Pick, curved

Cawthorne Hooks

165mm long, black finish



Code	Description
Sb-ENT/206a	Slight curve
Sb-ENT/206b	Medium curve
Sb-ENT/206c	Full curve
Sb-ENT/206d	Double curve
Sb-ENT/206e	90° angle, 0.5mm
Sb-ENT/206f	90° angle, 1.0mm
Sb-ENT/206g	45° angle, 0.5mm
Sb-ENT/206h	Retrograde
Sb-ENT/206i	Straight

PICKS AND SCOOPS



Musgrove Cerumen Pick and Wool Carrier

Double ended



Code	Description
Sb-ENT/024	178mm long

Formby Cerumen Scoop and Hook

Double ended



Code	Description
Sb-ENT/026	178mm long

Hovell Loop


Fine, for removing hard wax



Code	Description
Sb-ENT/027	152mm long

Lempert Scoop

Oval, 200mm long



Code	Description
Sb-ENT/066a	2mm x 2.5mm
Sb-ENT/066b	2mm x 3.5mm
Sb-ENT/066c	2mm x 4mm
Sb-ENT/066d	2mm x 4.5mm
Sb-ENT/066e	3mm x 5.5mm
Sb-ENT/066f	3mm x 6mm

House Curette

Double ended, 180mm long



Code	Description
Sb-SH/5213	Straight
Sb-SH/5213a	Angled

Balance Scoop and Curette



Code	Description
Sb-ENT/038	230mm long

HOOKS AND PROBES



Hartmann Hook

Sharp, straight shaft



Code	Description
Sb-ENT/021	140mm long


Gardiner Brown Hook and Imray Loop

Double ended



Code	Description
Sb-ENT/025	127mm long

Cawthorne Cerumen Hook



Code	Description
Sb-ENT/023	143mm long

St Bart's Hook

Blunt, straight handle



Code	Description
Sb-ENT/022	178mm long

Tracheal Hook

Sharp



Code	Description
Sb-ENT/SH.5672	165mm long

McIndoe Skin Hook

Fine, with flat handle



Code	Description
Sb-ENT/057	191mm long

Tracheal Single Hook Retractor

Sharp



Code	Description
Sb-GI/SH5671	152mm long

Cawthorne Probe

Silver, flattened end



Code	Description
Sb-ENT/019	140mm long

Tracheal Double Hook

Blunt



Code	Description
Sb-ENT/SH.5673	105mm long

Jobson Horne Probe

Double ended, serrated tip and smooth ring




Code	Description
Sb-ENT/020a	140mm long
Sb-ENT/020b	178mm long

RETRACTORS

STAHLÉAN BRAZIL
SAVING LIVES


Mollison Retractor

Self retaining, curved



Code	Description
Sb-ENT/034a	4-4 teeth, adult, 165mm long
Sb-ENT/034b	


Wallace Modified Plester Retractor



Code	Description
Sb-ENT/062	2 x 2 teeth, 140mm long

Plester-Jansen Retractor


Solid blade one side,
140mm long



Code	Description
Sb-Ent/061A	Left
Sb-Ent/061B	Right
Sb-Ent/061C	3 x 3 Teeth

Shah Intra Aural Retractor

Satin finish
152mm long



Code	Description
Sb-ENT/060a	4-4 teeth, adult, 165mm long
Sb-ENT/060b	Third blade
Sb-ENT/060c	Intra aural blade
Sb-ENT/060d	Post aural blade

Hurd Pillar Retractor and Tonsil Dissector



Code	Description
Sb-ENT/147	197mm long

Mollison Anterior Pillar Retractor and Tonsil Enucleator

Semi sharp



Code	Description
Sb-ENT/148	210mm long

SUCTION TUBES



Ipswich Pattern Sucker



Code	Description
Sb-ENT/050a	267mm long

Lempert Suction and Irrigation Tubes


With tucker valve and stilette for cleaning
stainless steel, 203mm



Code	Description
Sb-ENT/044	6 F.G. 2mm
Sb-ENT/044a	9 F.G. 3mm

Formby Suction Tube

Angled, conical, tip size 10 F.G.
stainless steel



Code	Description
Sb-ENT/042	210mm long

Magill Suction Tubes


Set of four with conical ends, one curved, three angled,
handle, stainless steel, 210mm



Code	Description
Sb-ENT/045	Various sizes
Sb-ENT/046a	3.0mm ø
Sb-ENT/046b	3.5mm ø
Sb-ENT/046c	4.0mm ø

Barron Suction Tube

Small size, 4 F.G.
Not for use with Zoellner ends



Code	Description
Sb-ENT/041	229mm long

Zoellner Suction Tube


With stilette, stainless steel



Code	Description
Sb-ENT/047	170mm long, 2mm ø


House Suction and Irrigation Tubes

Stainless steel, 152mm long



Code	Description
Sb-ENT/043	8 F.G. x 6 F.G.
Sb-ENT/043a	10 F.G. x 7 F.G.
Sb-ENT/043b	12 F.G. x 8 F.G.

Rowbotham Suction Tube



Code	Description
Sb-ENT/169	286mm long

SPECULUM


STAHLÉAN BRAZIL
SAVING LIVES



Rosen Slotted Speculum

Set of four, round ends, black finish
38mm long


Code	Description
Sb-ENT/049j	Set of 4
Sb-ENT/049f	4mm
Sb-ENT/049g	5mm
Sb-ENT/049h	6mm
Sb-ENT/049i	7mm



Tumarkin Speculum

Set of eight, tapered slot, black finish
50mm long


Code	Description
Sb-ENT/052	Set of 8
Sb-ENT/053a	4.5mm ø
Sb-ENT/053b	6.5mm ø
Sb-ENT/053c	7.0mm ø
Sb-ENT/053d	7.5mm ø
Sb-ENT/053e	8.0mm ø
Sb-ENT/053f	8.5mm ø
Sb-ENT/053g	9.0mm ø
Sb-ENT/053h	11.0mm ø



Rosen Aural Speculum

Round, 40mm long

Code	Description
Sb-ENT/049a	3mm
Sb-ENT/049b	4mm
Sb-ENT/049c	5mm
Sb-ENT/049d	6mm
Sb-ENT/049e	7mm



Heath Speculum

Set of four, black finish
50mm long

Code	Description
Sb-ENT/207e	Set of 4
Sb-ENT/207a	4.0mm ø
Sb-ENT/207b	5.0mm ø
Sb-ENT/207c	6.0mm ø
Sb-ENT/207d	7.0mm ø




Shea Speculum

Set of five, oblique ends, black finish, 38mm long

Code	Description
Sb-ENT/050	Set of 5
Sb-ENT/051	3mm
Sb-ENT/051a	4mm
Sb-ENT/051b	5mm
Sb-ENT/051c	6mm
Sb-ENT/051d	7mm



Peters Seigle Speculum



Code	Description
Sb-ENT/366	Peters ear speculum black chrome finish set of 5, sizes 0 to 4
Sb-ENT/189i	2.2mm x 3.2mm
Sb-ENT/189j	2.8mm x 4.2mm
Sb-ENT/189k	4.3mm x 5.8mm
Sb-ENT/189l	5.2mm x 7.6mm
Sb-ENT/189m	6.0mm x 7.6mm
Sb-ENT/189d	6mm
Sb-ENT/189e	7mm
Sb-ENT/189f	8mm
Sb-ENT/189g	9mm
Sb-ENT/189h	10mm


SPECULUM

STAHLÉAN BRAZIL
SAVING LIVES



St Clair Thompson Nasal Speculum


Code	Description
Sb-ENT/072	38mm long
Sb-ENT/072a	51mm long
Sb-ENT/072b	64mm long
Sb-ENT/072c	89mm long



Thudichum Nasal Speculum

Set of seven, sizes 1-7,
70mm long

Code	Description
Sb-ENT/073a	Standard
Sb-ENT/073b	Insulated



Standard

Code	Description
Sb-ENT/073a	Size 1 - 12mm blade
Sb-ENT/073a	Size 2 - 14mm blade
Sb-ENT/073a	Size 3 - 16mm blade
Sb-ENT/073a	Size 4 - 18mm blade
Sb-ENT/073a	Size 5 - 20mm blade
Sb-ENT/073a	Size 6 - 22mm blade
Sb-ENT/073a	Size 7 - 24mm blade
Sb-ENT/073a	Size 8 - 27mm blade


Insulated

Code	Description
Sb-ENT/073b	Size 1 - 12mm blade
Sb-ENT/073b	Size 2 - 14mm blade
Sb-ENT/073b	Size 3 - 16mm blade
Sb-ENT/073b	Size 4 - 18mm blade
Sb-ENT/073b	Size 5 - 20mm blade
Sb-ENT/073b	Size 6 - 22mm blade
Sb-ENT/073b	Size 7 - 24mm blade
Sb-ENT/073b	Size 8 - 27mm blade

Killian Nasal Speculum

With Stop


Code	Description
Sb-ENT/081a	51mm blade
Sb-ENT/081b	64mm blade
Sb-ENT/081c	76mm blade
Sb-ENT/081d	89mm blade



Without Stop


Code	Description
Sb-ENT/082b	76mm blade

Lempert Endaural Speculum



Code	Description
Sb-ENT/060e	133mm long


Pilcher Nasal Speculum



Code	Description
Sb-ENT/071	152mm

Cottle Nasal Speculum

127mm long



Code	Description
Sb-ENT/192c	20mm blade
Sb-ENT/192b	32mm blade
Sb-ENT/192	41mm blade
Sb-ENT/192a	60mm blade


AURAL FORCEPS



Triangular Jaw



Smooth Jaw



Serrated Jaw



Tenaculum Jaw



Cup Jaw



Grunwald Punch



Henckel Jaw

Beale Aural Forceps

Micro, black finish, triangular flat jaws, serrated



Code	Description
Sb-ENT/183	70mm to shoulder

Hough Aural Forceps

Micro, extra fine oval cup jaws, black finish



Code	Description
Sb-ENT/059	70mm to shoulder

Cawthorne Aural Forceps

Micro, black finish, extra fine
70mm to shoulder



Code	Description
Sb-ENT/182	Smooth jaws
Sb-ENT/199	Serrated jaws

Ormerod Aural Forceps

Micro, extra fine, black finish
70mm to shoulder



Code	Description
Sb-ENT/184a	Serrated, angled upwards
Sb-ENT/184b	1mm cup jaw
Sb-ENT/184c	Elongated cup jaw, 0.5 x 4mm

Hartmann Aural Forceps

Micro, black finish
75mm to shoulder



Code	Description
Sb-ENT/058	Serrated jaws
Sb-ENT/058a	Tenaculum jaws
Sb-ENT/058b	Cup shape jaws
Sb-ENT/058c	Curved left, serrated
Sb-ENT/058d	Curved right, serrated
Sb-ENT/058e	Angled up, serrated
Sb-ENT/058f	Serrated jaws, satin finish
Sb-ENT/058s	Pre-sterile - single use only

Shea Aural Forceps

Micro, oval cup jaw, black finish
70mm to shoulder




Code	Description
Sb-ENT/069	Standard
Sb-ENT/069a	Curved to left
Sb-ENT/069b	Curved to right
Sb-ENT/069c	Curved upward

AURAL FORCEPS



Hartmann Aural Forceps


Standard, satin finish
75mm to shoulder



Code	Description
Sb-ENT/181	Serrated jaws
Sb-ENT/181a	Tenaculum jaws
Sb-ENT/181b	Cup shape jaws
Sb-ENT/181c	Grunwald punch action

Heath Aural Forceps


Light, with lifting discs, serrated jaws



Code	Description
Sb-ENT/175	83mm to shoulder

Henckel Aural Forceps


Fine fenestrated scoop blades



Code	Description
Sb-ENT/176	75mm to shoulder

Struycken Cutting Forceps

75mm to shoulder




Code	Description
Sb-ENT/178	Fenestrated under jaw

ARTERY FORCEPS



Birkett Tonsil Artery Forceps



Fine, 178mm long



Code	Description
Sb-ENT/151a	Straight
Sb-ENT/151b	Curved


Negus Tonsil Artery Forceps

180mm long



Code	Description
Sb-ENT/197	Large curve
Sb-ENT/197a	Small curve


Denis Browne Tonsil Artery Forceps



Code	Description
Sb-ENT/155	216mm long

Sawtell Tonsil Artery Forceps


178mm long



Code	Description
Sb-ENT/152	Large curve
Sb-ENT/153	Small curve

Fickling Artery Forceps


190mm long



Code	Description
Sb-PL/044b	1-2 teeth
Sb-PL/044a	Serrated

Scott Tonsil Artery Forceps

Curved to side, screw joint




Code	Description
Sb-ENT/154	180mm long

TONSIL, TURBINATE AND SEPTUM FORCEPS



Denis Browne Tonsil Holding Forceps


Angled shanks, no rack



Code	Description
Sb-ENT/144	200mm long

Jansen-Middleton Nasal Septum Forceps


Compound action
178mm long



Code	Description
Sb-PL/032	5mm jaw


Irwin Moore Nasal Turbinate Forceps

191mm long




Code	Description
Sb-ENT/077d	5mm jaw

Walsham Septum Forceps



Code	Description
Sb-PL/030a	Straight, 215mm long
Sb-PL/030b	Right, 215mm long
Sb-PL/030c	Left, 215mm long

Luc Nasal Turbinate Forceps



Heart shaped, 178mm long


Code	Description
Sb-ENT/078a	Small 6mm ø
Sb-ENT/078b	Medium 7.5mm ø
Sb-ENT/078c	Large 9.0mm ø

Oval, 191mm long

Code	Description
Sb-ENT/078d	Small 6mm ø
Sb-ENT/078e	Medium 7.5mm ø
Sb-ENT/078f	Large 9.0mm ø

Hartmann Conchotomes

Circular cutting
102mm to shoulder




Code	Description
Sb-ENT/080a	5mm ø
Sb-ENT/080b	7mm ø
Sb-ENT/080c	9mm ø
Sb-ENT/080d	11mm ø

DRESSING FORCEPS



Hartmann Dressing Forceps


Serrated jaw, screw joint



Code	Description
Sb-ENT/076	108mm to shoulder

Wilde Dressing Forceps


Angular, serrated jaws



Code	Description
Sb-ENT/075	102mm to shoulder

Tilley Dressing Forceps


Straight, fine, serrated ends



Code	Description
Sb-ENT/077	2mm wide, 76mm to shoulder

Irwin Moore Thyroid Clamp

Thin jaws, angled to side




Code	Description
Sb-ENT/170	165mm long

OTHER FORCEPS




Trousseau Tracheal Dilator



Code	Description
Sb-ENT/172	Adult 140mm long
Sb-ENT/172a	Child 125mm long

Wilson Bone Nibbling Forceps


Fine, curved, satin finish
152mm long



Code	Description
Sb-ENT/054	Downward cutting - 3mm jaws
Sb-ENT/055	Upward cutting - 3mm jaws


Cawthorne Nibbling Forceps

Narrow, 3mm jaws, wide handle



Code	Description
Sb-ENT/068	159mm long

Denis Browne Collar Crushing Forceps



Code	Description
Sb-4062	140mm long

Lempert Nibbling Forceps

152mm long




Code	Description
Sb-ENT/067a	Straight, 2.5mm jaws
Sb-ENT/067b	Angled, 3.0mm jaws

MOUTH GAGS

STAHLÉAN BRAZIL
SAVING LIVES




Doyen Mouth Gag



Code	Description
Sb-ENT/88d	Adult
Sb-ENT/88e	Child


Fergusson Ackland Mouth Gag



Code	Description
Sb-ENT/209	Adult, 210mm
Sb-ENT/209a	Child, 155mm


Boyle-Davis Mouth Gag

Frame only



Code	Description
Sb-ENT/083a	Adult
Sb-ENT/083b	Child


Mason Mouth Gag



Code	Description
Sb-ENT/210	Adult, 210mm


Kilner Dott Mouth Gag

Frame only



Code	Description
Sb-ENT/087a	Adult
Sb-ENT/088e	Set, including frame, small, medium, and large blades

Draffin Bipod



Adult

Code	Description
Sb-ENT/083	Adult, 4 hole x 480mm long half circle end for use with a drape
Sb-ENT/083E	Adult, 4 hole x 480mm long with round ends to be used with Magauran plate
Sb-ENT/083E	Adult, 4 hole x 480mm long with round ends to be used with Magauran plate

Child

Code	Description
Sb-ENT/083C	Child, 3 hole x 425mm long half circle end for use with a drape
Sb-ENT/083D	Child, 3 hole x 425mm long with round ends to be used with Magauran plate
Sb-ENT/083F	Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern
Sb-ENT/083D	Child, 3 hole x 425mm long with round ends to be used with Magauran plate
Sb-ENT/083F	Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern
Sb-5104	Magauran plate 480mm - for use with Draffin bipods


Magauran Plate

Code	Description
Sb-5104	480mm - for use with Draffin bipods

TONGUE PLATES AND DEPRESSORS



Doughty Tongue Plate



Plain

Code	Description
Sb-ENT/084a	65mm - 1.6mm-1.9mm slot
Sb-ENT/084b	75mm - 1.6mm-1.9mm slot
Sb-ENT/084c	90mm - 4mm-4.5mm slot
Sb-ENT/084d	
Sb-ENT/084e	

Slotted


Code	Description
Sb-ENT/085a	65mm - 1.6mm-1.9mm slot
Sb-ENT/085b	75mm - 1.6mm-1.9mm slot
Sb-ENT/085c	90mm - 4mm-4.5mm slot
Sb-ENT/085d	
Sb-ENT/085e	115mm - 4mm-4.5mm slot

With anaesthetic tube

Code	Description
Sb-ENT/086a	65mm - 1.6mm-1.9mm slot
Sb-ENT/086b	75mm - 1.6mm-1.9mm slot
Sb-ENT/086c	90mm - 4mm-4.5mm slot
Sb-ENT/086d	100mm - 4mm-4.5mm slot
Sb-ENT/086e	115mm - 4mm-4.5mm slot

Kilner Dott Tongue Plate

With endotracheal tube groove



Code	Description
Sb-ENT/088a	Small
Sb-ENT/088b	Medium
Sb-ENT/088c	Large
Sb-ENT/088d	Extra small


Tongue Depressor

Flat stainless steel



Code	Description
Sb-ENT/143	184mm long

Lack Tongue Depressor




Code	Description
Sb-ENT/142a	Child, 13mm wide
Sb-ENT/142b	Adolescent, 19mm wide
Sb-ENT/142c	Adult, 22mm wide

PUNCHES




Nasal Cutting Forceps



Myles

Through cutting, 120mm to shoulder


Code	Description
Sb-ENT/SH5889	12.5 x 3.5mm jaw



Gruenwald-Henke

Through cutting, 120mm to shoulder


Code	Description
Sb-ENT/SH5281	7 x 2mm jaw
Sb-ENT/SH5282	9 x 2.5mm jaw



Struycken


Through cutting, 120mm to shoulder

Code	Description
Sb-ENT/SH5280	11 x 1mm jaw



Grunwald Nasal Punch Forceps


108mm to shoulder



Code	Description
Sb-ENT/077a	A fenestrated jaw, 3.5mm
Sb-ENT/077b	Fenestrated lower jaw, narrow, 2.5mm
Sb-ENT/077c	Fenestrated and tapered lower jaw, 2-3mm

Seymour Jones Antrum Punch Forceps


159mm long



Code	Description
Sb-ENT/132	4mm jaw

Backbiter


180mm long



Code	Description
Sb/4765b	7 x 3mm jaw

Backbiter

200mm long




Code	Description
Sb-ENT/5392	7 x 2mm jaw

PUNCHES



Stammberger Rhinoforcep Antrium Punch


7mm x 3mm, 130mm long



Code	Description
Sb-ENT/234	Left back cut
Sb-ENT/234a	Right back cut

Tilley Henckel Ethmoid Punch


Graduated shank, 114mm to shoulder



Code	Description
Sb-ENT/138a	4.5mm
Sb-ENT/138b	5.5mm
Sb-ENT/138c	6.5mm
Sb-ENT/138d	8mm

Ostrom Antrum Punch


Forcep, 178mm to shoulder



Code	Description
Sb-ENT/135	3mm jaw

Mackay Modified Stammberger Punch

200mm long



Code	Description
Sb-ENT/5392	7 x 2mm jaw

Paterson Biopsy Forceps

Elongated cup, 3mm jaw

Code	Description
Sb-ENT/174a	305mm long
Sb-ENT/174b	355mm long
Sb-ENT/174c	406mm long
Sb-ENT/174d	508mm long




PUNCHES




Chevalier Jackson Endoscope Forcep

With crocodile action and serrated jaw



Code	Description
Sb-ENT/241	Effective shaft length 300mm
Sb-ENT/242	Effective shaft length 400mm
Sb-ENT/243	Effective shaft length 500mm
Sb-ENT/244	Effective shaft length 600mm




Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and curved cup



Code	Description
Sb-ENT/245	Effective shaft length 300mm
Sb-ENT/246	Effective shaft length 400mm
Sb-ENT/247	Effective shaft length 500mm
Sb-ENT/248	Effective shaft length 600mm



Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and straight cup



Code	Description
Sb-ENT/249	Effective shaft length 300mm
Sb-ENT/250	Effective shaft length 400mm
Sb-ENT/251	Effective shaft length 500mm
Sb-ENT/252	




St Bartholomew's Ethmoid Punch

Scoop blades, 3mm bite, 114mm to shoulder



Code	Description
Sb-ENT/140a	Straight
Sb-ENT/140b	Upward cutting
Sb-ENT/140c	Downward cutting




PUNCHES

STAHLÉAN BRAZIL
SAVING LIVES



Ferris Smith Kerrison Punch


Heavy, 5mm bite,
229mm long



Code	Description
Sb-ENT/137a	Upward cutting
Sb-ENT/137B	Downward cutting

Citelli Antrum Punch


76mm to shoulder



Code	Description
Sb-ENT/133a	3mm jaw
Sb-ENT/133b	5mm jaw

Hajek Antrum Punch


197mm long



Code	Description
Sb-ENT/134	5mm jaw

Hajek Sphenoidal Punch

146mm to shoulder



Upward cutting

Code	Description
Sb-ENT/141a	4mm bite
Sb-ENT/141b	5mm bite

Downward cutting


Code	Description
Sb-ENT/141c	4mm bite
Sb-ENT/141d	5mm bite

RONGEURS



Beck Rongeur


Straight, for intervertebral disc



Code	Description
Sb-NEU/125	2mm bite, 130mm shaft

Pennybacker Rongeur


120mm shaft



Code	Description
Sb-NEU/133	3mm bite
Sb-NEU/107	4mm bite

Cushing Rongeur

130mm shaft




Code	Description
Sb-NEU/127-2	2mm bite
Sb-NEU/127-3	3mm bite
Sb-NEU/127-4	4mm bite
Sb-NEU/127-5	5mm bite
Sb-NEU/127-6	6mm bite



St Bart's Rongeur


114mm shaft



Code	Description
Sb-ENT/140a	3mm bite

Cushing Rongeur

130mm shaft



150° upward cutting


Code	Description
Sb-NEU/128c	2mm bite
Sb-NEU/128	3mm bite
Sb-NEU/128a	4mm bite
Sb-NEU/128b	5mm bite

150° downward cutting

Code	Description
Sb-NEU/130c	2mm bite
Sb-NEU/130	3mm bite
Sb-NEU/130a	4mm bite
Sb-NEU/130b	5mm bite

St Bart's Rongeur

140°, 3mm, 114mm shaft




Code	Description
Sb-ENT/140b	Upward cutting
Sb-ENT/140ac	Downward cutting

RONGEURS



Mackay Rongeur


140°, straight, 114mm shaft



Code	Description
Sb-SH/4249	5.5mm

Weil-Blakesley Rongeur

110mm shaft



140° upward cutting


Code	Description
Sb-ENT/194b	3.6mm bite
Sb-ENT/194e	4.2mm bite
Sb-ENT/194f	4.8mm bite

140° downward cutting

Code	Description
Sb-ENT/194g	3.6mm bite
Sb-ENT/194h	4.2mm bite
Sb-ENT/194i	4.8mm bite

Mackay Rongeur


With suction tube



Code	Description
Sb-SH/4249a	3.5mm bite, 114mm shaft

Weil-Blakesley Rongeur


120mm shaft



Code	Description
Sb-ENT/193b	3.0mm bite
Sb-ENT/193	3.6mm bite
Sb-ENT/194	4.2mm bite
Sb-ENT/194a	4.8mm bite
Sb-ENT/194c	5.8mm bite
Sb-ENT/194d	6.6mm bite

Mackay Rongeur

With suction tube, 140° upward cutting



Code	Description
Sb-SH/4249b	3.5mm bite, 114mm shaft

SCISSORS



Beckmann Turbinectomy Scissors


Curved



Code	Description
Sb-ENT/097	165mm long


Fomon Scissors

Straight, serrated blades, angled handles



Code	Description
Sb-ENT/200	145mm long

Scott Stevenson Tonsil Scissors



Code	Description
Sb-ENT/159	216mm long

Cottle Scissors


Straight blades, angled handles, blunt/blunt



Code	Description
Sb-ENT/OK201	165mm long

Wood Tonsil Scissors


Sharpened outer edges



Code	Description
Sb-ENT/156	178mm long

Cottle Bulldog Scissors


Blunt/blunt, 115mm long



Code	Description
Sb-ENT/202	115mm long

Wheeler Tonsil Scissors


Triangular form blades, serrated shanks



Code	Description
Sb-ENT/157	184mm long

Joseph Scissors


Sharp/sharp, 150mm long



Code	Description
Sb-ENT/203	Straight
Sb-ENT/204	Curved

Tilley Tonsil Scissors


Curved



Code	Description
Sb-ENT/158	190mm long

Heymann Turbinectomy Scissors

Angular




Code	Description
Sb-ENT/098	197mm long

SCISSORS



Wilson Tonsil Scissors



Code	Description
Sb-ENT/160	203mm long


Cawthorne Scissors

Micro, extra fine, 70mm to shoulder



Code	Description
Sb-ENT/186c	Left
Sb-ENT/186d	Right
Sb-ENT/186e	Straight

Thyroid Scissors



Code	Description
Sb-ENT/171	152mm long

Ormorod Aural Scissors


Micro, extra fine, curved upwards



Code	Description
Sb-ENT/185	70mm to shoulder

Knapp Scissors


Slight curve, blunt points, 114mm long



Code	Description
Sb-ENT/187a	Satin finish
Sb-ENT/187b	Black finish

Portmann Aural Scissors

Micro, satin finish 70mm to shoulder




Code	Description
Sb-ENT/186a	Curved to right
Sb-ENT/186b	Curved to left

GOUGES, KNIVES AND ELEVATORS



Ballenger Nasal Cartilage Knife

Revolving blade



Straight, 203mm long

Code	Description
Sb-ENT/096a	6mm wide
Sb-ENT/096b	8mm wide

Cranked, 197mm long

Code	Description
Sb-ENT/096c	6mm wide

Sickle Knife



Code	Description
Sb-BUPNS002	188mm long
Sb-ENT/271	160mm long

Jenkins Gouge

140mm long



Code	Description
Sb-ENT/037	1mm
Sb-ENT/037a	2mm
Sb-ENT/037b	3mm
Sb-ENT/037c	4mm
Sb-ENT/037d	5mm
Sb-ENT/037e	6mm
Sb-ENT/037f	7mm

Code	Description
Sb-ENT/037g	8mm
Sb-ENT/037h	9mm
Sb-ENT/037i	10mm
Sb-ENT/037j	11mm
Sb-ENT/037k	12mm
Sb-ENT/037l	13mm
Sb-ENT/037m	14mm

Killian Nasal Gouge

127mm long



Code	Description
Sb-ENT/093	4mm

Guy's Hospital Nasal Gouge

Shallow groove and flattened back, 127mm long



Code	Description
Sb-ENT/094	4mm

Tilley Nasal Gouge

V shaped, 127mm long



Code	Description
Sb-ENT/095	4mm

Cottle Elevator

Double ended, black chrome finish
200mm long



Code
Sb-ENT/224

Hill Elevator

Double ended, curved to side, rugine end



Code	Description
Sb-ENT/091	5mm wide, 222mm long

Tumarkin Atticotomy Elevator

Double ended




Code	Description
Sb-ENT/036	191mm long

CURETTES AND MYRINGOTOMES



Buck Ear Curette


Angled, blunt, 146mm long



Code	Description
Sb-ENT/028	2mm
Sb-ENT/028a	2.5mm
Sb-ENT/028b	3mm
Sb-ENT/028c	3.5mm
Sb-ENT/028d	4mm

Nasal Curette


Straight, 185mm long



Code	Description
Sb-ENT/628702	5 x 2mm

Beckmann Adenoid Curette


216mm long



Code	Description
Sb-ENT/167a	8mm blade
Sb-ENT/167b	10mm blade
Sb-ENT/167c	12mm blade
Sb-ENT/167d	14mm blade
Sb-ENT/167e	16mm blade
Sb-ENT/167f	18mm blade

Agnew (Politzer) Myringotome

190mm long



Code	Description
Sb-ENT/031	Shaft length 100mm


Trautmann Myringotome

178mm long



St Clair Thomson Adenoid Curette

With cage, 235mm long



Code	Description
Sb-ENT/168a	8mm
Sb-ENT/168b	10mm
Sb-ENT/168c	12mm
Sb-ENT/168d	14mm
Sb-ENT/168e	16mm
Sb-ENT/168f	18mm


Code	Description
Sb-ENT/032	Upward cutting
Sb-ENT/033	Downward cutting

ELEVATORS AND DISSECTORS



Freer Instruments

Single ended, 190mm long



Knife

Code	Description
Sb-ENT/065a	Curved to side
Sb-ENT/065c	Rounded edge, 8mm blade
Sb-ENT/065d	Flat rounded blade
Sb-ENT/065e	Rounded blade, curved on flat

Elevator

Code	Description
Sb-ENT/065f	Dull edged
Sb-ENT/065h	Sharp edged

Chisel

Code
Sb-ENT/065i

St Clair Thompson Dissector

159mm



Code	Description
Sb-ENT/089a	Sharp
Sb-ENT/089b	Blunt

Negus Aspiring Dissector


With 6mm olive, 210mm long



Code	Description
Sb-ENT/090a	-
Sb-ENT/090b	Insulated

Freer Separator and Elevator

Double ended



Code	Description
Sb-ENT/065b	Blunt/sharp, 180mm long
Sb-ENT/065g	Sharp, knife curved to side, 180mm long
Sb-ENT/065j	Blunt/sharp, 200mm long
Sb-ENT/065j-sal	Blunt/sharp, 188mm long

Syme Tonsil Dissector

Double ended, one end blunt, one semi sharp



Code	Description
Sb-ENT/149	170mm long

Gwynne-Evans Tonsil Dissector

Double ended, with serrated end and blunt end
197mm long




Code	Description
Sb-ENT/146a	6mm wide
Sb-ENT/146b	8mm wide

MISCELLANEOUS



Woakes Eustachian Catheters


Silver, Luer, 140mm



Code	Description	Code	Description
Sb-ENT/015	7 F.G.	Sb-ENT/015c	10 F.G.
Sb-ENT/015a	8 F.G.	Sb-ENT/015d	11 F.G.
Sb-ENT/015b	9 F.G.		

Kramer Eustachian Catheters


Silver, Luer, 140mm



Code	Description	Code	Description
Sb-ENT/014	7 F.G.	Sb-ENT/014c	10 F.G.
Sb-ENT/014a	8 F.G.	Sb-ENT/014d	11 F.G.
Sb-ENT/014b	9 F.G.		

Myles Naso-Antral Perforator


171mm long



Code	Description
Sb-ENT/103a	16 F.G.
Sb-ENT/103b	19 F.G.

Tilley Antrum Burrs

171mm long



Code	Description
Sb-ENT/131a	6mm ø
Sb-ENT/131b	7mm ø
Sb-ENT/131c	8mm ø

Laryngeal Mirror


Mirror only, 130mm long



Code	Description	Code	Description
Sb-ENT/191-6	6mm ø	Sb-ENT/191-20	20mm ø
Sb-ENT/191-8	8mm ø	Sb-ENT/191-22	22mm ø
Sb-ENT/191-10	10mm ø	Sb-ENT/191-24	24mm ø
Sb-ENT/191-12	12mm ø	Sb-ENT/191-26	26mm ø
Sb-ENT/191-14	14mm ø	Sb-ENT/191-28	28mm ø
Sb-ENT/191-16	16mm ø	Sb-ENT/191-30	30mm ø
Sb-ENT/191-18	18mm ø		

Laryngeal Mirror

Complete with metal handle, 260mm long



Code	Description	Code	Description
Sb-ENT/192-6	6mm ø	Sb-ENT/192-20	20mm ø
Sb-ENT/192-8	8mm ø	Sb-ENT/192-22	22mm ø
Sb-ENT/192-10	10mm ø	Sb-ENT/192-24	24mm ø
Sb-ENT/192-12	12mm ø	Sb-ENT/192-26	26mm ø
Sb-ENT/192-14	14mm ø	Sb-ENT/192-28	28mm ø
Sb-ENT/192-16	16mm ø	Sb-ENT/192-30	30mm ø
Sb-ENT/192-18	18mm ø		

Mitchell Trimmer



Code	Description
Sb-ENT/013	165mm

MISCELLANEOUS



Lake Seeker

Double ended, for vestibule



Code	Description
Sb-ENT/040	184mm long


Negus Pusher

"Knot tier", metal handle




Code	Description
Sb-ENT/165	203mm long

Lempert Rugine



Code	Description
Sb-ENT/064a	Narrow, 156mm long
Sb-ENT/065a	Broad, 171mm long

Stanford Cade Insufflator



Code	Description
Sb-ENT/188	175mm long

Irwin-Moore Tonsil Needle

203mm long



Code	Description
Sb-ENT/164a	Curved to right
Sb-ENT/164a	Curved to left

Tubing Mount

For use with catheters, Luer



Code	Description
Sb-ENT/016	25mm



STAHLÉAN BRAZIL

SAVING LIVES

CONTACT US

Office in são paulo

☎ +5511 972111007

☎ +92-3347711007

📷 @stahleanbrazil

✉ info@stahleanbrazil.com

🌐 www.stahleanbrazil.com